
EN	Dual Snap Lead Wire
FR	Fil de connexion enfichable double
DA	Strømførende ledning med dobbelt trykknop
NL	Aansluitkabel met dubbele knipsluiting
DE	Doppel-Elektrodenkabel mit Clip
IT	Filo di derivazione a doppio attacco
NO	Dobbel avledning som kan kneppes på
PT	Fio condutor de encaixe duplo
ES	Alambre guía con presión dual

86-41314-0010 Rev G.

EN

CAUTION STATEMENTS

Any serious incident that has occurred in relation to this device should be reported to Dymedix Diagnost and the competent authority of the Member States in which the user and/or patient is established.

Symbol	Definition of Symbol
	Latex Free
	Do not re-use
	Federal Law (USA) restricts this device to sale by or on the order of a physician.
	Medical Device
	Consult instructions for use.
	Date of Manufacture Expressed as: YYYYMMDD
	Manufactured by

FEATURES AND INTENDED USE:

The Dymedix Dual Snap Lead Wire has been designed to connect diagnostic snap electrodes to notice.

TECHNICAL SPECIFICATIONS

Connectors:

- Red and White matched connectors
- 1.5mm DIN socket safety connectors to the recording equipment
- 0.159" (4.03mm) dome snap connectors to the electrodes

Lead Wire:

- Two conductor, grey.

Resistance:

- End-to-end is less than 10 ohms

PRECAUTIONS FOR EFM CABLE

1. Store at room temperature.
2. Do not sterilize.
3. Connect to electrically isolated equipment.
4. Ensure wire and connectors are completely dry before using.
5. When removing the lead wire from electrodes or equipment, hold by the molded end and pull – DO NOT PULL BY THE WIRE. If wires or connectors become damaged, discontinue use.

CLEANING

1. Use a cleaning agent that is not corrosive to plastics.
2. Wipe the wire and connectors.
3. Do not immerse in liquids.

WARRANTY

Dymedix warrants its sensors to be free from defects in materials and workmanship and to perform satisfactorily for One (1) Year from the date of purchase. Should the electrodes fail, the sole liability of Dymedix is limited to replacement of the product with no charge for parts or labor. Under no circumstances shall Dymedix be liable for any consequential damages arising from the use of, or the inability to use this product. This warranty is rendered void for conditions resulting from: damage, marginal performance or malfunctions caused by: reuse, misuse, abuse, neglect, improper amplification settings, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original purchaser.

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Dymedix Diagnostics Inc.