Natus Disposable Adhesive (Surface) Electrodes

Instructions for Use:

Natus Neurology Incorporated
3150 Pleasant View Road
Middleton, Wisconsin 53562 USA
Tel: +1 608-829-8500
Website: natus.com

EU Representative / Importer
Natus Manufacturing Limited
IDA Business Park
Gort, Co. Galway, Ireland

Natus Nicolet UK Ltd
Baynards Green Trading Estate
Prospect House
Oxfordshire
Bicester
OX27 7SG
England, United Kingdom

Arazy Group Switzerland GmbH
Bruderholzallee 53
4059 Basel
Switzerland
swiss.ar@arazygroup.com

Rx Only

Associated product part numbers:

019-400400, 019-400500, 019-406600, 019-409000, 019-409100, 019-414200, 019-415000,
019-415200, 019-420800, 019-422200, 019-439300, 019-439400, 019-435300, 019-435500,
019-768500, N2025, N6003, 443010, 454010, 443020, 454030 and 454040

017W066 Rev 13 12/2022

English

Description:

Natus disposable adhesive (surface) electrodes are used as a non-invasive method of conducting electrical activity produced by muscles. Surface electrodes are used in stimulation and recording of electrical activity for analysis using an EMG amplifier.

Intended Use:

The intended use of surface recording electrodes is to stimulate or record the electrical activity to or from muscles and motor and sensory nerves during EMG nerve conduction studies, evoked potential studies, or inter-operative monitoring (IOM).

Intended User and Target Patient Group:

The product is for use only under the direction of trained medical professionals in acquiring electrophysiologic signals with patients selected at their discretion.

Clinical Benefits:

Facilitates obtaining an EMG recording to detect any irregularities indicative of muscle disorders.

Contraindications and Side Effects:

There are no known contraindications or side effects for procedures performed with Natus Disposable Adhesive (Surface) Electrodes. See Warnings or Precautions below.

Operating Instructions:

- Open the pouch just before use.
- Prepare the electrode site prior to application making sure that the site is clean and dry.
- Remove the electrode from its backing.
- Apply to the patient making sure to place even pressure across the entire electrode area.
- Re-seal pouch after each use.
- Electrodes can be used to record or stimulate.
- Following use, dispose of adhesive electrodes with infectious waste.

Understanding Cautions Statements:

 CAUTION
Refers to a hazardous situation that could result in minor or moderate injury or material damage if not avoided.
• Information on how the hazardous situation is avoided.

Warnings and Precautions:

 CAUTION
Device dropped or damaged in transit/use could lead to loss of function or delayed diagnosis. • Inspect the device prior to each use and do not use if damaged.
Electrode reused on another patient leads to cross infection or loss of performance. • Do not reuse or sterilize.
Data cannot be recorded due to improper placement leading to delay in procedure. • Assure proper connections and check signal quality before using.
Adhesive electrodes left in place for prolonged periods of time may be associated with skin redness. • Do not use the electrodes for prolonged periods of time.
Repositioning of electrode can degrade the signal quality. • Electrode should be replaced if the adhesive no longer provides a uniform connection to the skin surface.
Modifications to the device can affect function and performance. • Do not modify this equipment without authorization of the manufacturer.

Environmental Specifications:

Operating Conditions:

- Temperature: +10°C (+50°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 70 kPa to 106 kPa

Storage Conditions:

- Temperature: 0°C (+32°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 50 kPa to 106 kPa

Disposal Instructions:

Following use, dispose of adhesive electrodes with infectious waste.

Disclaimer:

Natus Medical Incorporated is not responsible for injury, infection or other damage resulting from the use of this product.

Any serious incident that has occurred in relation to the device should be reported to Natus Medical Incorporated and the competent authority of the Member State in which the user and/or patient is established.

Refer to the Natus website for an electronic copy of this document.

Glossary of Symbols

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
	-	-	An indication of Medical device	This product is a medical device.
Rx Only	21 CFR Part 801.109(b)(1)	Labeling-Prescription devices.	Prescription only	Indicates the product is authorized for sale by or on the order of a licensed healthcare practitioner.
	ISO 15223-1 Symbol 5.4.5 (Reference Annex B for the general prohibition symbol)	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Not made with Natural Rubber Latex	Indicates a medical device that is not made with natural rubber latex.
	-	-	Quantity	Number of parts in a package.

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
	ISO 15223-1 Symbol 5.1.1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Manufacturer	Indicates the medical device manufacturer.
	ISO 15223-1 Symbol 5.1.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Authorized representative in the European Community	Indicates the Authorized representative in the European Community.
	ISO 15223-1 Symbol 5.1.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Symbol 5.1.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Batch or Lot code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Symbol 5.1.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Symbol 5.4.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	ISO 60601-1 Table D.1 #11	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Operating instructions	
	ISO 15223-1 Symbol 5.4.3 Annex A #A.15	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Consult instructions for use	Indicates an instruction to consult an electronic instructions for use (eIFU).
	ISO 60601-1 Table D.2 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Follow instructions for use	Refer to instruction manual/ Booklet. NOTE on ME EQUIPMENT "Follow instructions for use"
	ISO 15223-1 Symbol 5.4.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Caution: Read all warnings and precautions in instructions for use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 60601-1 Table D.1 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.		
	ISO 60601-1 Table D.2 #2	Medical electrical equipment — Part 1: General requirements for basic safety and Essential performance.	General warning sign	Indicates a hazard of potential personal injury to patient or operator.
	ISO 15223-1 Symbol 5.1.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Use-by date	Indicates the date after which the medical device is not to be used. Note: This symbol shall be accompanied by a date to indicate that the medical device should not be used after the end of the year, month or day shown.
	ISO 15223-1 Symbol 5.2.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1 Symbol 5.3.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Temperature limit	Indicates the (storage) temperature limits to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.3.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Humidity limitation	Indicates the range of (storage) humidity to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.4.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Swiss Medical Device Ordinance (MedDO)	Swiss Medical Device (MedDO).	Indicates the Authorized Representative in Switzerland	Indicates the Authorized Representative in Switzerland.
	UKCA Medical Device Regulation (SI 2022 No 618, as amended) (UK MDR 2002)	UKCA Medical Device Regulation.	UKCA Mark	In compliance with the United Kingdom technical conformity.
	EU Medical Device Regulations 2017/745	EU Medical Device Regulation.	CE marking	Signifies European technical conformity.

Product Part Numbers	Label Part Numbers	Product Description
019-400400	031134	1.0m (39") lead (1.5mm DIN), 4/pouch
	031135	1.0m (39") lead (1.5mm DIN), 60 (15 sets of 4/box)

Product Part Numbers	Label Part Numbers	Product Description
019-400500	031130	1.0m (39") lead, (1.5mm DIN), 1/pouch
	031131	1.0m (39") lead, (1.5mm DIN), 24/box
019-406600	031136	1.0m (39") lead (1.5mm DIN), 2/pouch
	031137	1.0m (39") lead (1.5mm DIN), 48 (24 pairs/box)
019-409000	031138	2.0m (79") lead (1.5mm DIN), 4/pouch
	031139	2.0m (79") lead (1.5mm DIN), 60 (15 sets of 4/box)
019-409100	031132	2.0m (79") lead (1.5mm DIN), 1/pouch
	031133	2.0m (79") lead (1.5mm DIN), 24/box
019-414200	031140	1.0m (39"), (1.5mm DIN), 3/pouch
	031141	1.0m (39"), (1.5mm DIN), 72 (24 sets of 3/box)
019-415000	031142	1.5m (59") lead (1.5mm DIN), 4/pouch (2 red, 2 black)
	031143	1.5m (59") lead (1.5mm DIN), 96 (24 sets of 4/box)
019-415200	031144	1.0m (39") lead (1.5mm DIN), (1 Ground, 2 Disks), 3/pouch
	031146	1.0m (39") lead (1.5mm DIN), (1 Ground, 2 Disks), 72 (24 sets/box)
019-420800	031147	3.0m (118") leads (1.5mm DIN), 4/pouch (2 red, 2 black)
	031148	3.0m (118") leads (1.5mm DIN), (2 each red, black), 96 (24 sets of 4/box)
019-422200	031123	1.5m (59") lead, (1.5mm DIN), 1/pouch
	031124	1.5m (59") lead, (1.5mm DIN), 24/box
019-435300	031153	10/pouch Use with alligator-clip lead wires
	031154	100 (10 sets of 10/box) Use with alligator-clip lead wires
019-435500	031155	10/pouch Use with alligator-clip lead wires
	031156	100 (10 sets of 10/box) Use with alligator-clip lead wires
019-439300	031149	0.5m (20") lead, (1.5mm DIN), 2/pouch
	031150	0.5m (20") lead, (1.5mm DIN), 20 (10 sets of 2/box)
019-439400	031151	0.5m (20") lead, (1.5mm DIN), (1 each Red, Black, Green), 3/pouch
	031152	0.5m (20") lead, (1.5mm DIN), (1 each Red, Black, Green), 72 (24 sets of 3/box)
019-768500	031157	2/pouch Use with alligator-clip lead wires
	031158	50 (25 sets of 2/box) Use with alligator-clip lead wires
443010	031125	15mm x 20mm, 16.5cm (6.5") lead, (0.7mm male touchproof), 12/pouch
443020	031127	15mm X 20mm, 1.2m (47") lead, (1.5mm DIN), 12/pouch
454010	031126	20mm X 25mm, 16.5cm (6.5") lead, (0.7mm male touchproof), 12/pouch
454030	031128	20mm X 25mm, 1.2m (47") lead, (1.5mm DIN), 12/pouch
454040	031129	20mm X 25mm, 2.0m (79") lead, (1.5mm DIN), 12/pouch
N2025	031159	60 each Bar, Ring, and Ground Electrodes Use with alligator-clip lead wires
N6003	031160	12mm x 18mm contact area, 10/pouch Use with alligator-clip lead wires
	031161	12mm x 18mm contact area 100 (10 sets of 10/box) Use with alligator-clip lead wires

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