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Cortex strip and grid electrode

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Kortex Streifen-und Gridelektrode

Electrodo cortical

Eletrodo cortical

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English

☐ **Overview and use**

The cortex strip and grid electrode is an intracranial multiple-contact electrode used in functional neurosurgery for the electrophysiological exploration of the surface of the cortex.

The cortical electrode may be used in the following cases (see paragraph “Specifications”):

- presurgical diagnosis of certain forms of drug resistant epilepsy (period of use generally one week, and in any case limited to 30 days) or
- intraoperative mapping of functional areas for brain surgery (period of use limited to the duration of the surgery, i.e. less than 24 hours).

Use of the cortex electrode is strictly limited to these cases and must only be carried out in a specialised hospital setting.

☐ **Specifications**

Two cortex electrode configurations are available:

- Cortex electrodes references C10-nnXM (where nn is the number of contacts and X represents the number of contacts per strip)
These electrodes, comprising between 4 and 64 platinum/iridium contacts, with a usable diameter of 2.5 mm, are designed for use in presurgical diagnosis of certain forms of drug resistant epilepsy or for intraoperative mapping for brain surgery. They have connectors allowing sheath tunnelling for chronic implantation. They are shipped with blue connection holders and must be connected to recording / stimulation equipment via connection systems (see help sheet on “Cortical electrode connection”). The electrical resistance between a contact and the electrode connector is about 45 Ω; the insulation resistance between two contacts is more than 500 MΩ under 500 V.
- Cortex electrodes references C10-nnXIOM (where nn is the number of contacts and X represents the number of contacts per strip)
These electrodes, comprising between 4 and 16 stainless steel contacts, with a usable diameter of 4 mm (shouldered contacts) may be used only for intraoperative mapping for brain surgery. They have a connection allowing them to be connected directly to recording / stimulation equipment (cable 1.8 m long with female plugs diameter 1.5 mm to comply with standard DIN 42802). The electrical resistance between a contact and the electrode connector is around 500 Ω; the insulation resistance between two contacts is more than 500 MΩ under 500 V.

Cortex electrodes are available in the form of “strips” or “grids”. The grids are pre-cut into strips of 4, 6 or 8 contacts, thus the neurosurgeon is completely free to choose the configuration best suited to the patient. The electrode contacts are numbered to facilitate intraoperative exploration.

The materials in contact with the patient are silicone and, depending on the specific cortex electrode used, either platinum/iridium or stainless steel.

☐ **Symbols guide**

	“Catalogue reference”
	“Batch number”
	“Single-use item”
	“Sterilised using an ethylene oxide process”
	“Use by”
	“Important”
	“Serial number ”
	“Refer to Instructions”
	“Manufacturer”
	“Store in a dry place”
	“Do not resterilize”
	“Keep away from sunlight”
	“Do not use if the packaging is damaged”

☐ **Warnings and recommendations**

⇒ **Recommendations prior to use**

The cortex electrode is shipped in sterile condition for single use. Before use, it is necessary to check the integrity of the packaging and device.

⇒ **Positioning of the electrode**

The implantation of a cortex electrode requires a specific neurosurgical procedure which can only be carried out under the supervision of an experienced neurosurgeon (see paragraph on *Typical implantation procedure*).

Warning: Due to their fragile nature, cortex electrodes must be handled with care. In particular, **no excessive pulling force should be applied to the electrode sheaths**.

Case of cortex electrodes C10-nnXM:

A flexible spatula may be needed to position the electrode beyond the bone flap. During chronic implantation, the electrode sheaths must be tunnelised to minimise the risk of infection.

Case of cortex electrodes C10-nnXIOM:

These electrodes are designed to be placed on the surface of the cortex, in the craniotomy zone, during brai surgeries. **They should not be positioned beyond the bone flap**, between the bone and the brain.

Warning: Cortex electrodes C10-nnXIOM must not be used for chronic implantation.

⇒ **Electrode connection**


The electrode may only be connected to systems specifically suited to electroencephalographic (EEG) exploration, and in compliance with currently applicable safety requirements.

Case of cortex electrodes C10-nnXM:

DIXI Microtechniques offers connexion systems with a coding device to allow quick identification of each contact, thus avoiding any risk of confusion.

The DIXI Medical cortex electrode has been developed for connection with a DIXI Medical connexion system, and the DIXI Medical connexion systems have been developed for connection to a DIXI Medical electrode. DIXI Microtechniques declines all responsibility in the case of a DIXI Medical electrode being used with a non-DIXI Medical connexion system, or a DIXI Medical connexion system being used with a non-DIXI Medical electrode.

To enable the electrode to be connected, each electrode connector must be held in a blue connection holder shipped with the electrode (the end of the connector must be lined up with one end of the blue connection holder).

Warning: Once in position, the blue connection holder must not be removed from the electrode connector, as this may damage the connection.

Warning: The connection holders shipped with the electrodes are essential for their connection. Make sure they are transferred to the recording room when they are not positioned on the electrode connectors directly after implantation. In case of problems, DIXI Microtechniques offers sets of 10 blue connection holders SCC-08S-10.

An electrode with its connection holder is then connected:

- either to a DIXI Medical extension cable ref. PCP-16* or PCU-16*, which is itself connected to a DIXI Medical adaptor ref. ADP-*,
- or to a DIXI Medical connection cable ref. CBL-1x* or CBL-2x*.

For more information about these connection systems, please refer to their respective Instructions for use.

Case of cortex electrodes C10-nnXIOM:

These electrodes are designed for direct connection to recording / stimulation equipment. The cortex electrode connector plugs are numbered to allow quick identification of each contact, thus avoiding any risk of confusion.

Warning: The electrode connection is fragile. At the end of surgery, **the electrode connector plugs must be disconnected by pulling on the plug**. Pulling on the sheath risks breaking off the connector plug.

⇒ **Use with other devices**

To date, there has been no contra-indication concerning the use of other devices normally found in a neurosurgery block environment (defibrillator, electrosurgical devices etc.) at the same time as the cortex electrodes. However, the practitioner is responsible for initiating suitable monitoring procedures where specific devices are to be used.

⇒ **Exposure to Magnetic Resonance Imaging (MRI)**

Case of cortex electrodes C10-nnXM:

The platinum / iridium alloy used for cortex electrode conductors is a non ferromagnetic material with low magnetic susceptibility.

However, the compatibility of cortex electrodes with an MRI examination has not been demonstrated. The subjection of a patient with implanted electrodes to an MRI scan is the responsibility of the neuroradiologist, upon condition that the following precautions are taken:

- 1/ avoid the formation of loops in the cable to minimise the risk of inducted current.
- 2/ for the same reasons, make sure that there can be no contact made between the different implanted electrode connectors or with other conductive parts. The blue DIXI Medical connection holder must be positioned on the electrode connector
- 3/ when interpreting the data, take into account the risk of distortion or disruption of the images
- 4/ do not connect the electrodes to any device

Warning: In all cases, this examination must be carried out under the supervision of a neuroradiologist and limited to emergency situations or when postoperative location of the electrodes is necessary.

Case of cortex electrodes C10-nnXIOM:

These cortical electrodes contain materials which are not compatible with MRI examination. DIXI Microtechniques declines all responsibility concerning the conduct of MRI scans when these electrodes are used.

⇒ **Contra-indications**

Cortex electrodes must not be used in the following cases:

- General contra-indications for any surgical operation
- Infectious cutaneous lesions to the scalp
- Anticoagulant, platelet antiaggregant or non-steroidal anti-inflammatory treatment

⇒ **Side effects**

The surgical procedure and the use of cortex electrodes can lead to the following complications (complications issued of scientific literature):

- Infectious complications such as meningitis, cerebral abscess, epidural abscess, subdural empyema, subdural empyema, osteomyelitis, scalp abscess, scalp infection, subgaleal infection, cerebritis;
- Hemorrhagic complications such as subdural hematoma, epidural hematoma, intraparenchymale hemorrhage;
- Cerebrospinal fluid leak;
- Cerebral edema and mass effect of the electrodes;
- Increase of the intra-cranial pressure which can be linked with a supratentorial herniation, a pneumocephalus or a brain edema ;
- Neuro-deficitis such as hemiparesis, aphasia, sensory deficit, visual deficit, nerve palsy, Gerstmann syndrome, cavernous sinus syndrome, headache ;
- Removal or displacement of implanted electrodes ;
- Other medical complications that may not be direct consequences of cortex electrodes implantation: fever, headache, nausea, vomiting, need of blood transfusion (especially in the pediatric patients), anemia, muscle fibrosis, atrial fibrillation, deep venous thrombosis, pulmonary edema, pulmonary embolism, pneumonia, respiratory distress syndrome, acute renal failure, otitis, conjunctivitis, eye discharge, pilonidal cyst, cellulitis, diarrhea, urinary tract infection.

⇒ **Stimulation parameters**

To avoid any risk associated with heating of the electrode, the stimulation current must not exceed the limit value of 20 mA, at a maximum frequency of 200 Hz.

⇒ **Monitoring**

Continuous monitoring of the patient whose electrodes have been connected up to the EEG exploration systems is essential to avoid any risk of incidents.

⇒ **Disposal**

The cortex electrode is not recyclable and must be disposed of in accordance with standard contaminated hospital waste procedures.

☐ **Packaging contents**

Case of cortex electrodes C10-nnXM:

The sterile packaging contains a cortex strip or grid electrode and as many individually-wrapped blue connection holders as the electrode has connectors. An instruction for use and connection help sheet are enclosed in each box.

Case of cortex electrodes C10-nnXIOM:

The sterile packaging contains a cortex strip or grid electrode. An instruction for use is enclosed in each box.

☐ **Sterilisation and storage**

This device is shipped STERILE for SINGLE USE and must be used before the “EXPIRY DATE” indicated on the packaging. The cortex electrode has been sterilised with ethylene oxide.

DIXI Microtechniques cannot guarantee the levels of contamination after use, nor the functionality of the product in the case of re-use. **The cortex electrode should not be resterilised or re-used.**

The cortex electrode must be stored in its original packaging at room temperature, in a dry place, protected from physical impacts and light. It is absolutely essential that the packaging be checked before use.

Warning: If the seal is broken, do not use the device.

☐ **Typical implantation procedure**

Each neurosurgeon, in accordance with his own experience, will decide upon the most suitable implantation procedure. In general terms, the main stages are:

- Positioning of the electrodes on the surface of the cortex after opening the neurocranium and the meninges in accordance with the adopted surgery plan.

Case of cortex electrodes C10-nnXM when electrodes are used for chronic exploration

- Identification of the electrodes using the colour code on the electrode sheaths.
- Suture of the dura mater around the electrode sheaths.
- Tunnelling of the sheaths under the skin to reduce the risk of infection, using a needle (internal diameter > 1.65 mm, i.e. 13 Ga minimum).
- Closure of the bone flap
- Fixing of sheaths by sutures to the skin

Warning: the quality of the sutures around the electrode sheaths is very important in limiting cerebrospinal fluid loss.

☐ **Limited replacement and liability guarantee**

DIXI Microtechniques guarantees that the electrodes conform to requirements when used within the expiry date shown on the packaging. Items may not be returned without the prior approval of DIXI Microtechniques (returned electrodes are to be sent for inspection in appropriate packaging and include a detailed description of the defect). The warranty expressly excludes faults or defects resulting in part or in full from failure to comply with the storage conditions or from improper use.

No other warranty is given, whether expressly or implicitly, beyond the limited replacement warranty. DIXI Microtechniques declines all responsibility with regard to any other undertaking offered in its name in respect of the product and prohibits anyone from offering any other such undertaking.

DIXI Microtechniques declines all responsibility for any direct or consequential damage relating to the use of this product such as:

- external causes not linked to products supplied by DIXI Microtechniques or causes which are beyond its control.
- use of the product in a manner that does not strictly comply with the recommendations, instructions or warnings given previously herein.

Therefore, all consequence related to the reuse of the product, may not be attributed to DIXI Microtechniques (see paragraph *Sterilisation and storage*).

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