Cuffie per EEG - Istruzioni per l'uso.....	1
EEG Headcap - Instruction for use.....	6
EEG Elektrodenhaube - Gebrauchsanweisung.....	11
Bonnets EEG - Instructions.....	16
EEG Headcap - Brugsanvisning.....	21
Gorro EEG - Instrucción de utilizo.....	26
Touca para EEG - Instrução de uso	31
EEG kape - Upute za korištenje.....	36
Czepek do EEG - Instrukcja użytkowania.....	41
EEG Čepice - Návod k obsluze.....	46
脑电帽 - 使用说明书.....	51

Cuffie per EEG

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Headcap for EEG recordings.

EN ISO 13485

Manufacturer: Spes Medica S.p.A., Via Buccari 21, 16153 Genoa, Italy

All rights are reserved for this user manual, in particular the right to reproduction and distribution as well as translation. No part of this user manual may be reproduced in any form (e.g. photocopy, microfilm or other procedure) or processed, copied or distributed using electronic system without the prior written consent of Spes Medica.

CONTENTS

1. Intended Use

2. Handling

2.1 Packaging of the EEG Headcap

2.2 Instructions

3. Cleaning and disinfection

3.1 Cleaning

3.2 Deep disinfection

4. Waste disposal

5. Warnings

6. Symbols

7. Technical specifications

7.1 Material

7.2 Connectors and cables

8. Manufacturer

1. Intended use

Spes Medica EEG Headcaps are passive medical devices, used for EEG recordings. They have been developed to ensure maximum signal quality and patient comfort thanks to the combination of different materials and configurations that make them suitable for different customer needs.

Spes Medica EEG headcaps are made in accordance with the 10-20 International System up to 128-channels application with the Reference electrode (REF) integrated directly on the headcap.

2. Handling

The recorded signal is closely related to the right placement on the patient's head and an appropriate skin preparation.

2.1 Packaging of the EEG Headcap

EEG Headcap is provided in a plastic box, ready to be used.

2.2 Instructions for Use

Open the box and make a visual control to assure that the headcap is not damaged. Headcap positioning on the patient's head is a critical and essential step to get a good recorded signal.

After the correct placement of the headcap, it's necessary to secure it attaching a pair of elastic band to the chest band or using a chin holder.

Make sure that the headcap is stretched and adhered smooth.

The next step is to connect the headcap to the equipment.

Then is necessary to prepare the patient's skin carefully, otherwise the recorded signal may be affected by noise artifacts.

Therefore, it's necessary gently scrub the skin of the patient with a blunted needle to reduce skin impedance and finally inject enough amount of conductive gel to create a contact interface between the skin and the electrode.

In case of noisy signal or high impedance contact, verify again the previous steps.

3. Cleaning and disinfection

3.1 Cleaning

Spes Medica EEG Headcap needs to be cleaned immediately after every single use by carefully removing gel residues, water and mild detergent are required.

Once rinsed with water to remove all the residues of detergent, dry at room temperature without using other heat sources or other tools that could damage it.

3.2 Deep disinfection

For an appropriate disinfection, soak the headcap (non-critical use) in a solution, following the disinfectant data for concentration and time of action.

Spes Medica recommends the use of the commercial disinfectants shown in the table below:

Disinfectant	Active agent
Gigasept FF (new)	Succinate acid dialdehyde; dimethoxytetrahydrofuran
Mucocit- T	Didecyltrimethylammoniumchloride; alkylpropylenediamino-1,5-bisguanidiniumacetate; bis(-aminopropyl)-laurylamine; laurylpropylenediamine; non-ionic surfactants; auxiliary agents.
Gigasept pearls	Sodium percarbonate; Tetraacetylenethylenediamine.
Korsolex plus	N-(3-aminopropyl)-N-dodecylpropane, Didecyltrimethylammonium chloride

Spes Medica headcaps have resisted to the disinfectants list above, i.e. this method doesn't change the functional features of the product.

4. Waste disposal

Spes Medica EEG Headcaps must be disposed of in accordance with regulations in force for medical devices.

5. Warnings

Take preventive measures to avoid the direct contact between the headcap connector and any kind of liquid

Let the headcap dry at room temperature without using external heat source and/or other tools that could damage it.

Pay attention cleaning each headcap electrode by removing all traces of gel.

Make sure that the patient's skin under the headcap is not damaged.

6. Symbols and explanations

Symbols	Meaning
	Complies with the Regulation (EU) 2017/745
	Catalogue number
	Lot number
	Medical Device
	Consult instruction for use
	Warning
	Manufacturer
	Not Sterile

All Spes Medica EEG Headcaps with prewired electrodes have UDI (Unique Device Identification) barcode GS1-128.

7. Technical specifications

7.1 Materials

Spes Medica uses different combination of materials depending on the intent of use and the customer's specific requests.

In the table below the different materials used for the main components of the headcap are listed.

Component	Material
Headcap	Net, Neoprene

ElectrodeSn, Ag/AgCl, Ag/AgCl Sint, Au

Electrode holder

Silicon rubber, PE

7.2 Connectors and cables

Spes Medica offers different type of connectors depending on the user's equipment. The cables have also two possible configurations, flat cables and twisted cables, depending on the user's needs.

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