Instruction For Use

EN CNSAC Disposable EEG Cap

REF: EC-150, EC-151, EC-152, EC-163, EC-164, EC-165, EC-166, EC-167, EC-168, EC-169

Intended use: The Disposable EEG Cap is intended for use in routine clinical settings where rapid placement of a number of EEG electrodes is desired. Disposable EEG Caps are accessories for the diagnostic equipment (recording/monitoring device), and as such devices are intended to be used by qualified healthcare professional.

Related accessories: Disposable EEG Cap must be used in combination with CNSAC adapter cable and a DIN 42802 compatible EEG recording/monitoring device.

Application:
1.Clean the electrode position on the patient's scalp
2.Place the cap on the patient's scalp and fix it symmetrically on both sides of the head using the chin strap
3.Apply an appropriate amount of EEG conductive paste to the electrode and make sure that the electrode is placed completely on the scalp without shifting.
4.Connect the cap using cable adapter with the EEG recording device
5.Check teh signal quality of all electrodes before recording start
6.Start recording;
7.Once the recording is complete, remove the cap and help the patient to remove the paste from the skin
8.Dispose of the cap in accordance with hospital's regulations for medical waste

Product performance:
● Resistances≤100 Ω.
● AC impedances2 kΩ (less than 100 μA, peak-to- peak, at 10 Hz).
● DC offset voltages≤100 mV (after 1min stabilization period).
● Combined offset instability and internal noises≤150 μV
● Bias current tolerance (DC Voltage Offset) ≤100 mV

Contraindications:
● Patients with scalp or skin injuries in the area of application
● Patients who suffer from the negative psychological effects of a cap & hat

Cautions and warning:

1.The product should be used under the guidance of a healthcare professional.
2.Select the appropriate cap according to the size of the patient's head circumference and use EEG conductive paste and adapter cable to carry out the measurement
3.The product should be used immediately after opening the package
4.Check the product and packaging. If the product is damaged, please do not use it
5.Please do not use excessive force during use to prevent cables or electrodes from breaking
6.Do not use in the presence of flammable anesthetic gas mixed with air or flammable anesthetic gas mixed with oxygen or nitrous oxide
7.The disposable EEG cap is intended for single use and cannot be reused.

Disinfection and sterilization: The disposable EEG cap is non-sterile and intended for single patient use only. If the product needs to be disinfected, clean it with water, spray it with 75% alcohol (avoid connection) and then allow it to air dry. The use of strong disinfectant solutions and sterilization procedures may damage the product.

Storage conditions: Store in a temperature of -10°C ~+40°C, relative humidity ≤80%, atmospheric pressure 700hPa ~ 1060hPa. Free of corrosive gases, cool, dry, well-ventilated environment.

Note: Please report any serious incidents such as death, temporary or permanent serious deterioration of a patient, user, or third person, or a serious public health threat (in accordance with (EU) 2017/745 Annex 1, §23.4z) to the manufacturer and your national competent authority.

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