Genuine Grass® Gold Disc Electrodes

Instructions for Use:

Natus Manufacturing Limited
IDA Business Park
Gort, Co. Galway, Ireland
Tel.: +353 (0)91 647400
Website: natus.com

Arazy Group Switzerland GmbH
Bruderholzallee 53
4059 Basel
Switzerland
swiss.ar@arazygroup.com

UK Responsible Person
Optima Medical Ltd (Neuro).,
30a Alan Turing Road
The Surrey Research Park
Guildford Surrey
GU2 7AA United Kingdom

Rx only

Associated product part numbers:

FD-E5GH-48, FD-E5GH-60, FD-E5GH-72, FD-E5GH-96, F-E5GH-12, F-E5GH-120, F-E5GH-24, F-E5GH-30, F-E5GH-48,
F-E5GH-60, F-E5GH-96, F-E5GH-72, F-E6GH-48, FH-E5GH-48, FH-E5GH-60, FH-E5GH-72, FH-E5GH-96, FS-E5GH-120,
FS-E5GH-48, FS-E5GH-60, FS-E5GH-72 and FS-E5GH-96

English

Description:

Genuine Grass Gold Disc (cup) reusable electrodes are for surface monitoring EEG, electromyogenic and physiologic data signals commonly recorded on the chest, extremities, head and faces of patients. Cup electrodes are also used in electrocardiography, evoked potentials, and polysomnography. The Genuine Grass (R) Gold Disc Electrodes are MR (Magnetic Resonance) Unsafe.

Intended Use:

Surface recording of EEG, LTM, EP, EMG, ECG, EOG and other bioelectric signals.

Intended User and Target Patient Group:

Reusable Grass electrodes are for professional use only and should be used in compliance with accepted industry standards.

Clinical Benefits:

Facilitates obtaining the recording of an EEG/PSG study to detect any irregularities indicative of various brain or sleep disorders.

Contraindications and Side Effects:

There are no known contraindications or side effects for procedures performed with Genuine Grass Gold Disc electrodes.

Operating Instructions:

- Using a cotton tip applicator apply a small amount of skin prep gel such as NuPrep® or to the site of the electrode application and rub gently.
- Fill the electrode cup with a conductive paste or gel such as Ten20®. Do not overfill the cup. Excessive amount of paste or gel can result in asymmetric signals.
- Fix the cup to the patient with a small amount of pressure.
- Use a small pre-cut gauze square or piece of paper tape to help secure the cup electrode to the application site.
- Check that clear, strong signals are being transmitted.

Cleaning Instructions:

Cleaning:

Clean electrodes by immersion in warm, soapy tap water (15-25°C) using a mild detergent for a minimum of 1 minute. Do not immerse the touchproof connector in cleaning solution. Immersing the entire unit may reduce the life of the connector and should therefore be avoided. Rub electrodes gently using a soft brush to remove visible material. Dry with a clean, lint-free cotton cloth or air dry. Note that 24k gold plate is extremely soft and scratches easily.

Sterilization:

Note: Sterilization efficacy is the user's responsibility. The following methods are for reference only and Natus does not guarantee that effective sterilization will result.

Steam Sterilization:

Steam sterilization set points:

Sterilizer Type: Prevacum
Preconditioning Pulses: 4
Temperature: 132°C
Full Cycle Time: 4 minutes
Dry Time: 30 minutes
Configuration: Pouched

Operating Data	
Shelf life	1 year
Period of use cleaning	≤1040 cleanings
Period of use sterilization	≤200 sterilization cycles
Connector repeatability	200 mating/un-mating cycles

Sterilizer Type: Prevacum

Preconditioning Pulses: 4

Temperature: 134°C

Full Cycle Time: 3 minutes

Dry Time: 30 minutes

Configuration: Pouched

Ethylene Oxide (EO) Sterilization:

Condition for 30 minutes use standard settings. Exposure parameters:

Temperature: 55°C

Relative Humidity: 50%

EO Concentration: 725 mg/L

Gas Exposure Time (full cycle): 1 hour

Gas Type: 100% EO

Aeration time post-sterilization: 12 hours

STERRAD® Sterilization using the 100S and 100NX Systems with APTIMAX® Trays:

100S STERRAD System:

Place electrodes in an APTIMAX® tray and process in a STERRAD 100S sterilization half-cycle with the following parameters:

Pre-Plasma:	10 minutes
Injection:	6 minutes
Diffusion:	2 minutes
Post-Plasma:	2 minutes
Injection Volume:	1440 µL of 59% H ₂ O ₂

100NX STERRAD System:

Place electrodes in an APTIMAX tray and process in a STERRAD 100NX sterilization half-cycle with the following parameters:

Transfer:	8 minutes
Diffusion:	0.5 minutes
Plasma:	7.5 minutes
Injection Volume:	4900 µL of 53% H ₂ O ₂

Understanding Cautions Statements:

 CAUTION
Refers to a hazardous situation that could result in minor or moderate injury or material damage if not avoided.
• Information on how the hazardous situation is avoided.

Warnings and Precautions:

 CAUTION
Device dropped or damaged in transit/use could lead to loss of function or delayed diagnosis.
• Inspect the device prior to each use and do not use if damaged.
Device reused on another patient leads to cross infection or loss of performance.
• Do not reuse without sterilizing the device.

 CAUTION
Data cannot be recorded due to improper placement leading to delay in procedure. Assure proper connections and check signal quality before using.
Electrodes left in place for prolonged periods of time may be associated with skin redness. Do not use the electrodes for prolonged periods of time.
Modifications to the device can affect function and performance. Do not modify this equipment without authorization of the manufacturer.

Environmental Specifications:

Operating Conditions:

- Temperature: +10°C (+50°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 70 kPa to 106 kPa

Storage Conditions:

- Temperature: 0°C (+32°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 50 kPa to 106 kPa

Compliance Standards:

- ISO 10993: 2018 Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process
- ETS 300 019-2-1 Environmental Engineering (EE); Environmental conditions and environmental tests for Storage
- ASTM D4169-16 Standard Practice for Performance Testing of Shipping Containers and Systems for Vibration
- EN 60601-1:2006/A1:2013 and IEC 60601-1:2005/A1:2012 Ed. 3.1 clause 8.5.2.3

Disposal Instructions:

Natus is committed to meeting the requirements of the European Union WEEE (Waste Electrical and Electronic Equipment) Regulations 2014. These regulations state that electrical and electronic waste must be separately collected for the proper treatment and recovery to ensure that WEEE is reused or recycled safely. In line with that commitment Natus may pass along the obligation for take back and recycling to the end user, unless other arrangements have been made. Please contact us for details on the collection and recovery systems available to you in your region at natus.com.

Electrical and Electronic Equipment (EEE) contains materials, components and substances that may be hazardous and present a risk to human health and the environment when WEEE is not handled correctly. Therefore, end users also have a role to play in ensuring that WEEE is reused and recycled safely. Users of electrical and electronic equipment must not discard WEEE together with other wastes. Users must use the municipal collection schemes or the producer/importers take-back obligation or licensed waste carriers to reduce adverse environmental impacts in connection with disposal of waste electrical and electronic equipment and to increase opportunities for reuse, recycling and recovery of waste electrical and electronic equipment.

Equipment marked with the below crossed-out wheeled bin is electrical and electronic equipment. The crossed-out wheeled bin symbol indicates that waste electrical and electronic equipment should not be discarded together with unseparated waste but must be collected separately.

Disclaimer:

Natus Manufacturing Limited is not responsible for injury, infection or other damage resulting from the use of this product.

Any serious incident that has occurred in relation to the device should be reported to Natus Manufacturing Limited and the competent authority of the Member State in which the user and/or patient is established.

Refer to the Natus website for an electronic copy of this document.

Glossary of Symbols

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
	ISO 15223-1 Symbol 5.7.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	An indication of Medical device	This product is a medical device.
Rx only	21 CFR Part 801.109(b)(1)	Labeling-Prescription devices.	Prescription only	Indicates the product is authorized for sale by or on the order of a licensed healthcare practitioner.
	ISO 15223-1 Symbol 5.1.1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Manufacturer	Indicates the medical device manufacturer.

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
	ISO 15223-1 Symbol 5.1.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Symbol 5.1.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Batch or Lot code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Symbol 5.1.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Symbol 5.4.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	ISO 60601-1 Table 0.1 #11	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Operating instructions	
	ISO 60601-1 Table 0.2 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Follow instructions for use	Refer to instruction manual/ Booklet. NOTE on ME EQUIPMENT "Follow instructions for use"
	ISO 15223-1 Symbol 5.4.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Caution: Read all warnings and precautions in instructions for use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 60601-1 Table 0.1 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.		
	ISO 15223-1 Symbol 5.3.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Temperature limit	Indicates the (storage) temperature limits to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.3.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Humidity limitation	Indicates the range of (storage) humidity to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.2.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1 Symbol 5.4.5 (Reference Annex B for the general prohibition symbol)	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Not made with Natural Rubber Latex	Indicates a medical device that is not made with natural rubber latex.
	ISO 15223-1, Clause 5.2.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Non sterile	Indicates a medical device that has not been subjected to a sterilization process.
	2012/19/EU	Waste Electrical and Electronic Equipment (WEEE)	Disposal at end of operating life instructions	Indicates that electrical and electronic equipment waste should not be discarded together with unseparated waste but must be collected separately.

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
	ISO 15223-1 Symbol 5.1.2	Swiss Medical Device Ordinance (MedDO)	Indicates the Authorized Representative in Switzerland	Indicates the Authorized Representative in Switzerland.
	UK MDR 2002	UKCA Medical Device Regulation (SI 2002 No 618, as amended)	UKCA Mark	In compliance with the United Kingdom technical conformity.
	EU Medical Device Regulations 2017/745	EU Medical Device Regulation.	CE marking	Signifies European technical conformity.
	ASTM F2503	MR (Magnetic Resonance) Unsafe	MR Unsafe	Indicates that the product is not safe to be used in an MR environment.
	ISO 15223-1 Symbol 5.1.11	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Country of Origin	Indicates the country of origin.

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