Sliding Ruler

Règle à coulisse
Riga con cursore
Schiebelineal
Regla deslizante
Régua de deslizamento
Συρόμενος χάρακας
Плъзгаща се линия
Skydelære
Liukuva viivoitin
Räknesticka
Tolómérő
Skyvelær
Suwmiarka
Riglă glisantă
Posuvné měřítko
Bīdāmais lineāls
Posuvné meradlo
Лінійка з рухомим краєм

EN INSTRUCTION FOR USE	FI KÄYTTÖOHJE
FR MANUEL D'INSTRUCTION ET D'UTILISATION	SE BRUKSANVISNING
IT MANUALE DI ISTRUZIONI E PER L'USO	HU HASZNÁLATI UTASÍTÁS
DE ANLEITUNGS- UND BENUTZUNGSHANDBUCH	NO BRUKSANVISNING
ES MANUAL DE INSTRUCCIONES Y EMPLEO	PL INSTRUKCJA UŻYTKOWANIA
PT MANUAL DE INSTRUÇÃO E UTILIZAÇÃO	RO INSTRUCȚIUNI DE UTILIZARE
GR ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	CZ NÁVOD K POUŽITÍ
BG ИНСТРУКЦИИ ЗА УПОТРЕБА	LV LIETOŠANAS INSTRUKCIJA
DK BRUGSANVISNING	SK NÁVOD NA POUŽITIE
UA ІНСТРУКЦІЯ З ВИКОРИСТАННЯ	

2A route de Pouligny Chaudefontaine
25640 MARCHAUX-CHAUDFONTAINE France
Tel.: +33 (0) 3 81 88 98 90 - Fax: +33 (0) 3 81 88 98 99
info@diximedical.com

Date of revision: 04/2025
DIXI Medical template: EN-AR-16_R00
Date of obtaining CE marking: 1998

English

Intended use

The sliding ruler is used in stereotactic neurosurgery to set the depth of use of the drill bits and of the coagulation electrodes, and to measure the setting length of marking rods.

Its use is strictly reserved for this purpose and should only be carried out by specialized centers.

Medical indications

SEEG is indicated for evaluation of drug-resistant epilepsy, when a surgical hypothesis exists, but when anatomo-electroclinical data collected during the non-invasive phase are insufficiently concordant regarding the supposed localization of the region responsible for seizures (epileptogenic zone) and/or its relationship with functional areas. It is indicated to confirm the hypothesis and to establish surgical strategy.

The thermocoagulation function of the DIXI Medical depth electrode can also be used as a palliative treatment option in case of contraindication to standard treatment options for drug-resistant epilepsy (surgical resection).

Description

The sliding ruler:

- Is sold non-sterile for multiple use,
- Is autoclave-resistant up to 134° C,
- Is made of stainless steel,
- Enables to check the implantation depth of the accessories used (drill bits, coagulation electrodes, marking rod),
- Is graduated every millimeter,
- Has a measurement accuracy of ± 0,5 mm,
- Has a maximum measuring depth of 230 mm.

Clinical performances

The sliding ruler is an accessory used for the placement of DIXI Medical depth electrodes. It has an indirect clinical performance linked to DIXI Medical depth electrodes. The sliding ruler associated with other DIXI Medical accessories allows safe and efficient surgical placement of DIXI Medical Depth Electrodes, including subsequent electrodes removal, and connection to recording/stimulation/RF systems.

Clinical benefits

The sliding ruler is an accessory used for the placement of DIXI Medical depth electrodes. It has indirect benefits linked to DIXI Medical depth electrodes.

No direct clinical benefit, understood as a positive impact of the device on the health of individual, is expected from the use of the DIXI Medical depth electrodes alone. The indirect clinical benefit of the device usage lies through brain recording, stimulation and thermocoagulation capabilities, supporting more accurate medical decision-making by healthcare professionals. DIXI Medical depth electrodes can only be implanted and connected with the proprietary accessories. Given that accessories have no direct clinical benefit, but rather contribute to successful placement and connection of the electrodes, they therefore fall under the same defined indirect clinical benefit.

Symbol definition

	Medical device		Store in a dry place
	Manufacturer		Keep away from sunlight
	Consult instruction for use		Do not use if the packaging is damaged
	Caution		Temperature limit
	Catalogue number		Humidity limitation
	Serial number		Refer to instruction manual/booklet
	Date and country of manufacture		Federal (USA) law restricts this device to sale by or on the order of a physician.

Warnings and recommendations

Recommendations before use

The sliding ruler should be cleaned and sterilized prior to each use, even for the first use. The device being reusable, the hospital is responsible for the counting of the microbial load prior to sterilization, for the sterilization, for the checks for sterility after sterilization and for the sterile packaging. Do not use in case of blatant degradation of the device (for example a bad graduation readability) after cleaning and / or sterilization.

Target population

The device is intended for drug-resistant epilepsy adult patients and children from 1 year old and above 10 kg, who may be potential candidates for surgical resection. Pregnant women, breastfeeding women, neonates, premature neonates, and infants are excluded from the target population.

User profile

Level of training: The devices are used exclusively by medical teams in centers specialized in functional and stereotactic neurosurgery:
- Experienced neurosurgeons for the implantation, the use, and the explantation of the electrodes with the aid of their accessories.
- Hospital staff during the cleaning and sterilization process of DIXI Medical depth electrodes' reusable accessories.

Required knowledge:
SEEG is a procedure requiring specific equipment and a clinical team specialized in this technique. It should only be performed in hospitals with a functional and stereotactic neurosurgery department equipped to perform surgical treatment of epilepsy. Neurosurgeons are trained in using the Microdeep depth electrode and its accessories made by DIXI Medical.

Language comprehension:
Each device is provided with labelling and a user manual that complies with current regulatory requirements.

Experience:
The DIXI Medical depth electrodes and their accessories should be used by experienced and trained neurosurgeons, neurologists and neurophysiologists in the diagnosis and treatment of epilepsy.

Contraindications

The sliding ruler is without contact with the patient. There are no contraindications specific to the use of the sliding ruler to be declared. Nevertheless, the sliding ruler is used for the placement of DIXI Medical depth electrodes.

The DIXI Medical depth electrodes must not be used in the following cases:

- Presence of an infectious lesion of the skin;
- Anticoagulant and / or platelet antiaggregant treatment;
- General contraindications for any surgical operation such as fever, angina. .;
- Fine or fragile bones, which may limit the resistance of the guiding screw.

No thermocoagulation shall be performed near vascular structures (<2 mm from the geometric center of the selected contacts).

Side effects

The sliding ruler is without contact with the patient. There is no undesirable side-effect directly related to its use. Nevertheless, the sliding ruler is used for the placement of DIXI Medical depth electrodes.

The surgical procedure and the use of depth electrodes in an SEEG procedure may lead to the following complications (complications based on the scientific literature).

Common (≥1/100 to <1/20):
- Hemorrhagic complications such as hematoma (intracerebral, subdural, epidural), hemorrhage, bleeding or clots (intracerebral, subarachnoid, subdural...);
- Transient headache.

Uncommon (≥1/1000 to <1/100):
- Infectious complications such as intracranial infection, meningitis, abscess, scalp infection;
- Hydrocephalus;
- Cerebral edema;
- Pneumocephalus;
- Transient or permanent neurological deficit (aphasia, paresis, ...).

Rare (≥1/10000 to <1/1000):
- Cerebrospinal fluid leak;
- Encephalitis;
- Other medical complications: psychotic attack, status epilepticus, neurogenic pulmonary edema, arterial dissection, atelectasis, leg compartment syndrome, hyponatremia.

Not known (cannot be estimated from the available data):
- Other medical complications: subdural hygroma, cerebral infarction, cranial fracture, fever, leukocytosis, vasospasm, cardiac arrest, pneumonia, deep vein thrombosis, pulmonary embolism, urinary infection, gastroenteritis.

The thermocoagulation procedure may result in the following complications (complications based on the scientific literature):

Uncommon (≥1/1000 to <1/100):
- Transient or permanent neurological deficit (paresthesia, dysesthesia, hemianopsia, verbal memory deficit, apraxia, hemiparesis, ...);
- Slight sensation of heat or pain during coagulation.

Residual Risks

Failure to follow this instruction and use manual may result in burns, patient injury due to improper positioning of accessories, infectious complications, and prolonged procedure time.

Serious incident


Any serious incident that has occurred in connection to the device should be reported to the manufacturer and the competent authority in which the user and / or the patient is established.

Use of the device

The drill bit, the coagulation electrode and the marking rod must be handled with care when inserted into the sliding ruler. The depth of implantation must be determined beforehand by the neurosurgeon.

The sliding ruler is composed of two parts:

- A graduated ruler (1): it includes a groove enabling the insertion of the different devices.
- A slide (2) (sliding on the ruler): it is easily adjustable thanks to a screw that allows to move it along the ruler. The depth is measured at the edge of the slide. The slide fits perfectly the shape of the ruler. It is therefore inserted in the recess and acts as a stop for the devices.

- Interaction with other medical devices**

The sliding ruler should only be used with:

- The drill bits,
- coagulation electrodes and
- the marking rod.

The ACS-740 marking rod should be equipped with a stop to best adjust its depth of use. DIXI Medical offers stops adapted to the different accessories. It is recommended to use only DIXI Medical stops with medical devices made by DIXI Medical. For more information on stops please refer to the user manual.

DIXI Medical cannot be held responsible for any consequences arising from the use of sliding rulers with devices other than those it markets.

- Exposure to Magnetic Resonance Imaging (MRI)**

The sliding ruler is made of stainless steel. It should not be exposed to MRI. However, regarding its use, it is not intended to be used in an MRI environment.

- Device disposal**

This device is not recyclable and must be integrated to the standard treatment process of contaminated products in hospitals.

- ☐ **Packaging content**

The packaging contains a non-sterile sliding ruler packed in a bag. An instruction for use is enclosed in each packaging.

- ☐ **Cleaning and sterilization**

Warnings prior to the cleaning/sterilization procedure
The device is supplied non-sterile for multiple use. It should be cleaned and sterilized prior to each use according to the procedure described below.
The below-mentioned cleaning and sterilization procedure enable a full inactivation of non-conventional transmissible agents (NCTA) according to the Instruction DGS/R13/2011/449 of the 1st December 2011. Neurosurgical interventions which use DIXI Medical medical devices are acts posing a risk regarding NCTA. It is in any case necessary to consult the domestic instructions of the user country regarding the reduction of transmission risks of NCTA during invasive acts (Instruction DGS/R13/2011/449 of the 1st December 2011 for France).
This is the only procedure to be validated by DIXI Medical with cleaning/disinfection products and the specified parameters shown in italics in the instructions below. Moreover, the equipment used (washer-disinfector and autoclave) should be certified by the users and compliant with current standards (ISO 15883 for washer-disinfectors). Carrying out this procedure with other products or parameters of use, as well as the use of unqualified equipment is the responsibility of the user and must be carried out according to the <u>cleaning product manufacturers' instructions for use</u> .
The cleaning procedure, the microbiological contamination control (Bioburden), the sterilization method and the maintenance of the sterile state are the responsibility of the users.
PRE-DISINFECTION
Immersion
<u>Equipment</u> : no specific equipment. <u>Cleaning product</u> : enzymatic cleaning solution with the following properties: anti-corrosive (neutral pH, chloride-free), bactericidal, yeasticide, active on enveloped viruses. <u>Instructions</u> : Caution: carry out the points below within 30 minutes of using the device to prevent dirt present on the device from drying. - Prepare a cleaning product solution according to the supplier's instructions. - Completely immerse the device in the solution for the duration and at the temperature recommended in the supplier's instructions. <i>This step was validated by DIXI Medical using the product Aniosyme DD1 and aniosyme X3 (Laboratoires Anios) diluted to 0.5% for 5 minutes at room temperature.</i>
Rinsing
<u>Equipment</u> : no specific equipment. <u>Cleaning product</u> : water network. <u>Instructions</u> : Rinse the device in water network.
Drying
<u>Equipment</u> : clean material or a non-fuzzy gauze. <u>Cleaning product</u> : not applicable. <u>Instructions</u> : Dry the device up with a clean material or a non-fuzzy gauze. Following pre-disinfection, two solutions may be used: manual cleaning or automatic cleaning.
MANUAL CLEANING
Immersion
<u>Equipment</u> : no specific equipment. <u>Cleaning product</u> : cleaning agent (neutral pH) with the following properties: aldehyde free, anti-corrosive, bactericidal, yeasticide, mycobactericide, virucidal, active on prions and spores of bacteria. <u>Instructions</u> : - Prepare a cleaning product solution according to the supplier's instructions. - Clean the device with the solution prepared for the duration and at the temperature recommended in the supplier's instructions. - If necessary, eliminate residues in brushing the device gently with a soft plastic brush (Caution: The use of metal brushes is forbidden during this cleaning). <i>This step was validated by DIXI Medical using an undiluted Actanios HLD solution (Laboratoires Anios) for 30 minutes at room temperature (about 20°C).</i>
Rinsing
<u>Equipment</u> : no specific equipment. <u>Cleaning product</u> : sterile or filtered water (0,2μ). <u>Instructions</u> : Rinse the device with sterile or filtered water (0.2μ).
Drying
<u>Equipment</u> : clean material or a non-fuzzy gauze. <u>Cleaning product</u> : not applicable. <u>Instructions</u> : Dry the device up with a clean material or a non-fuzzy gauze.
AUTOMATIC CLEANING

Cleaning
<u>Equipment</u> : washer disinfector compliant to the ISO 15883 series. <u>Cleaning product</u> : alkaline cleaning agent (basic pH) with the following properties: anti-corrosive, bactericidal, yeasticide, mycobactericide, active on enveloped viruses, virucidal. <u>Instructions</u> : - Prepare a cleaning product solution according to the supplier's instructions. - Clean the device in a washer-disinfector with the solution prepared for the duration and at the temperature recommended in the supplier's instructions. <i>This step was validated by DIXI Medical using the Actanios LDI product (Laboratoires Anios) diluted to 1% for 10 minutes at 55°C.</i>
Rinsing
<u>Equipment</u> : washer disinfector compliant to the ISO 15883 series. <u>Rinsing product</u> : rinsing solution with acidic pH, which neutralizes alkaline residues of the cleaning agent and eliminates mineral residues from the water. <u>Instructions</u> : - Prepare a rinsing product solution according to the supplier's instructions. - Carry out a first rinse of the device in the washer-disinfector with the solution prepared for the duration and at the temperature recommended in the supplier's instructions. <i>This step was validated by DIXI Medical using the Anios RN product (Anios Laboratories) diluted to 0.2% for 2 minutes at 51°C.</i>
Rinsing
<u>Equipment</u> : washer disinfector compliant to the ISO 15883 series. <u>Rinsing product</u> : rinsing solution with acidic pH, drying activator. <u>Instructions</u> : - Prepare a rinsing product solution according to the supplier's instructions. - Carry out a second rinse of the device in the washer-disinfector with the solution prepared for the duration and at the temperature recommended in the supplier's instructions. <i>This step was validated by DIXI Medical using the Anios RHW product (Anios Laboratories) diluted to 0.1% for 3 minutes at 93°C.</i>
Drying
<u>Equipment</u> : washer disinfector compliant to the ISO 15883 series. <u>Cleaning product</u> : not applicable. <u>Instructions</u> : Dry the device in a washer-disinfector for a minimum of 15 minutes. <i>This step has been validated by DIXI Medical for 15 to 20 minutes at 105°C.</i>
PACKAGING
Packaging before sterilization
<u>Equipment</u> : no specific equipment. <u>Instructions</u> : place the device in packaging adapted to hospital practices (wraps, bags, etc.) and protecting it from shock.
STERILIZATION
Two steam sterilization cycles have been validated by DIXI Medical according to the ISO 17665 standard. The user can choose one of the following two cycles:
Cycle#1
<u>Time</u> : 18 minutes. <u>Temperature</u> : 134°C. <u>Sterilizer</u> : autoclave (steam)
Cycle#2
<u>Time</u> : 3 minutes. <u>Temperature</u> : 134°C. <u>Sterilizer</u> : autoclave (steam)

- ☐ **Service life**

The sliding ruler is a reusable device. It should be cleaned, decontaminated, and sterilized prior to each use (see paragraph "Cleaning and sterilization"). A visual inspection of the product should be carried out before and after each use.

Do not reuse the sliding ruler if its appearance (graduation readability...) or functioning is questionable or seems to be modified.

The sliding ruler has been designed to resist to a maximum of 50 uses. It is necessary to change the device after 50 uses.

- ☐ **Storage**

The sliding ruler should be stored protected from shocks, away from sunlight in a dry place. In particular, the sliding ruler should not be exposed to storage conditions outside the following temperature and humidity ranges:

- Temperature: [-25°C; 60°C].
- Relative humidity: [15%; 90%].

- ☐ **Procedure of use**

- Use of the sliding ruler with marking rod:**
Measure the distance "proximal end of the in guide - proximal end of the guiding screw" by inserting the marking rod into the instrument guide until it touches the guiding screw.
Secure this length on the marking rod with the stop and measure the length with the sliding ruler. Insert the marking rod into the groove of the sliding ruler until it rests on the slide, corresponding to the selected depth.

- Use of the sliding ruler with drill bits:**
Determine the desired length of the drill bit: once known, adjust the depth on the ruler with the slide.
Insert the drill bit in the groove of the ruler until it rests on the slide, corresponding to the chosen depth.
Advance the stop on the drill bit to the start of the rule and tighten it.

- Use of the sliding ruler with coagulation electrodes:**
Add 5 mm to the drill bit measurement on the sliding ruler.
Insert the coagulation electrode into the groove of the sliding ruler and transfer the length.
Advance the stop on the coagulation electrode and tighten it.

- ☐ **Limited guarantee for replacement and liability**

DIXI Medical guarantees the conformity of the sliding ruler within a limit of 50 uses. Each return must be approved beforehand by DIXI Medical (each device must be accompanied by a precise explanation about the defect found and sent back for assessment in an appropriate packaging). However, the defects and faults resulting, totally or partially, from storage alteration or mishandling are expressly excluded from the guarantee.

No other guarantee than the limited guarantee for replacement is expressed or implicit. DIXI Medical declines any responsibility with respect to any other commitment made in their name by anybody and forbids anybody to make such commitment.

DIXI Medical declines any responsibility for any direct or indirect damages related to the use of this product, such as those related to the following:

- Causes external to the products provided by DIXI Medical, or which are beyond their control.
- A use which is not strictly compliant with the previously described recommendations, instructions, or precautions.

Thus, in particular, any consequence related to a faulty sterilization or the maintenance of the device in non-sterile conditions must not, in any case, be imputable to DIXI Medical.

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