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Not all part numbers are available in all markets.

ENGLISH

Description:

The XactTrace® Single Use Respiratory Effort System (XactTrace belts) is a respiratory effort sensor used in the diagnosis of sleep disordered breathing. The XactTrace Respiratory Inductive Plethysmograph (RIP) sensor measures changes in inductance and converts it into a digital signal that gives both qualitative and quantitative data. This sensor is more sensitive and reliable than typical respiratory effort sensors, especially in cases of paradoxical breathing. This product is for use with Natus-compatible amplifiers.

Intended Use:

The XactTrace® Single Use Respiratory Effort Belt System is intended to measure respiratory effort to assist in the diagnosis of sleep disorders or sleep related respiratory disorders. The respiratory effort signals measured are processed to provide electrical signals suitable for connection to the inputs of physiological recording equipment. The intended environments are hospitals, institutions, sleep centers, sleep clinics. The XactTrace® Single Use Respiratory Effort Belt System is intended for diagnostics purposes only and is not intended to be used as an apnea monitor.

Indication for Use:

The XactTrace® Single Use Respiratory Effort Belt System is intended to diagnose sleep disorders or sleep related respiratory disorders by the acquisition of physiological respiratory effort during Sleep Studies (Polysomnography or PSG) and Long-Term Monitoring in epilepsy (LTM) and in patients of all ages.

- Sleep apnea or another sleep-related breathing disorder. In this condition, breathing repeatedly stops and starts during sleep.
- Periodic limb movement disorder. In this condition, the patient involuntarily flexes and extends the legs while sleeping. This condition is sometimes associated with restless legs syndrome.
- Narcolepsy. In this condition, the patient experiences overwhelming daytime drowsiness and sudden attacks of sleep.
- REM sleep behavior disorder. In this condition, the patient acts out dreams during sleep.
- Unusual behaviors during sleep. The provider may perform a test if the patient participates in unusual activities during sleep, such as walking, moving around a lot of rhythmic movements.
- Unexplained chronic insomnia. In this condition, the patient consistently has trouble falling asleep or staying asleep and the provider may recommend polysomnography.

Intended User:

A technologist or physician trained in administering PSG studies with specialized training and certification (RPSGT).

Target Patient Group:

The XactTrace® Respiratory Effort System can be used on patients of all ages. The size of the patient and duration of study determines the proper belt to use.

Clinical Benefit:

The XactTrace® Respiratory Effort Belt System measures respiratory effort to assist in the diagnosis of sleep disorders or sleep related respiratory disorders. The respiratory effort signals measured are processed to provide electrical signals suitable for the inputs of physiological recording equipment.

The XactTrace® Respiratory Effort Belt System is intended for diagnostic purposes only and is not intended to be used as an apnea monitor. Use the XactTrace® Respiratory Effort Belt System only under the direction and supervision of a physician or trained technologist.

Contraindications and Side Effects:

There are no contraindications identified for these supplies.

Operating Instructions:

Fitting the XactTrace® Belts

The XactTrace® belts must be custom fitted for each patient. Keep the XactTrace® belts dry during use.

1. Encircle the belt around the patient's chest under the arm to approximate the circumference for the thoracic belt.
2. When cutting the belt, reduce its circumference by 10-15cm (4-6 in) so it stretches around the thorax. The belt should fit snugly to prevent slippage during the study. It is important to use sharp scissors for a clean cut. The wire should not exceed the end of the belt.
3. Secure the cut ends of the belt into the belt lock with the blue connector according to the following steps:
 - Twist the top end of the belt lock clockwise to open the catch. The white mark on the top should match the open lock symbol on the catch.
 - Insert the cut end of the belt into the catch. Make sure to insert the end all the way to the bottom of the catch.
 - Twist the top end of the belt lock counter clockwise to close the catch. When the white mark on the top matches the closed lock symbol on the catch, the belt lock is properly closed.
4. Prepare the abdominal belt in the same manner, fitting the belt around the patient's stomach at the navel, and use the belt lock with the yellow connector.

Attaching the XactTrace® Belts to the Patient

The XactTrace® belts are intended to be worn over the patient's nightclothes.

1. Place the belt with the yellow connector around the patient's stomach, at the navel.
2. Take the two components of the belt lock and connect them.
3. Place the second belt with the blue connector around the patient's chest, under the arm.
4. Connect the second belt lock.
5. External inductive interface cable modules are required only with non-Natus amplifiers, referred to as "universal" applications. When using the inductive interface cables, insert the connector from the Chest belt into the processor module labeled THORAX, and then insert the connector from the Abdomen belt into the processor module labeled ABD. Both the Thorax (blue) and Abdomen (yellow) belt locks and inductive interface cable labels are color coded for easy identification and connection.
6. Plug the thorax and the abdomen module outputs into the appropriate bipolar inputs on the polygraph. Insert the red touch-proof plugs into (+) input, and the black plugs into (-) input.

User Settings

- **Sensitivity** – Sensitivity adjustment up or down is expected. Response depends on sensor application and patient effort.
- **Low Frequency Filter/Time Constant** – 0.03 Hz (or 5 seconds or longer)
- **High Frequency Filter** – 15 to 35 Hz.

Note: Shorter time constants or higher low frequency filter settings significantly attenuate waveforms.

Warnings and Cautions:

WARNING

A **WARNING** indicates that there is a risk of death or serious injury to the user or patient.

CAUTION

A **CAUTION** indicates that there is a risk of injury to the user or patient or risk of damage to the device.

WARNING

- Connect the XactTrace® belts only to an input that is electrically isolated from the mains power. Do not plug the cables into electrical outlets, as this could result in serious electrical shock.
- Take care not to cut any cables when cutting the belts.
- Do not use the device in an explosive environment, in the presence of flammable liquids, such as an anesthetic mixture with air, or with oxygen or nitrous oxide.
- Do not use damaged belts, sensors, or cables.

CAUTION

- Over-tightening the belt lock could sever the belt or result in impaired signal quality. Do not overtighten the belt locks.
- The belts are single use.
- Do not stretch the belts too tightly around the patient as this may cause discomfort.
- Care must be taken to ensure that cables do not encircle the patient's neck. Special attention is needed in the case of children.
- Device when used by untrained user could lead to delay in diagnosis. This device is intended to be used by qualified healthcare professionals.
- Avoid all unnecessary contact with moisture when using the device.
- This product is for diagnostic purposes only and is not to be used as an apnea monitor.
- Portable and mobile RF communications can affect the performance of the device.
- Electrostatic discharges (ESD) may cause artifacts in the signal from the device. A trained operator should be able to recognize these artifacts easily. Avoid conditions where electrostatic charge can build up due to low humidity and friction against carpets, clothing, and sheets made from artificial fibers.
- Do NOT use two abdomen or two thorax modules in the same recording. This can cause interference between the two sensors, which can result in unusable signals.
- Cables wrapped around the belt locks may cause the wire to break due to excessive strain on the cable. Do not wrap the cables around the belt locks.
- Cables wrapped around the sensor may cause the wire to break due to excessive strain on the cable. Do not wrap the cable around the sensor.
- XactTrace® belts shall be worn over the patient's clothing.
- Unauthorized modification could lead to loss of device function and performance. Do not modify device.

Environmental Specifications

	Operating conditions:	Storage and Transportation Conditions:
Temperature:	40°F to 120°F (+5°C to +50°C)	0°F to 120°F (-18°C to +50°C)
Relative Humidity:	15 to 95% (non-condensing)	10 to 95% (non-condensing)
Atmospheric Pressure	50 kPa to 106 kPa	50 kPa to 106 kPa

After completing use with a single patient, dispose of the belts and carefully store the XactTrace belt locks between recordings. Wrap the sensor cable loosely around the belt lock to protect the cable from damage. Disconnect modules from the belt locks when not in use to save battery life.

Output Specifications (Inductive Interface Cable)

Maximum signal amplitude:	< 2mV
Frequency Range:	0.2 to 3 Hz
Sensitivity:	Approximately 50 µV/mm

Cleaning and Disposal

Wipe the belt locks and inductive interface cable (if used) clean with a commercial wipe such as CaviWipes™ or Sani-Chlor® to remove visible soil. Wipe the article using a lint-free cloth and air dry. Take care to avoid contacting the sensor connector and plug of the belt locks or processor modules with the cleaning solution. The belt lock and inductive interface cable modules (if used) do not require sterilization. The belts should be discarded after use, and do not require cleaning.

Cable and Sensor Disposal

Cables and sensors can be used for one year from the date of purchase. Natus Medical Incorporated is committed to meeting the requirements of the European Union WEEE (Waste Electrical and Electronic Equipment) Directive 2012/19/EU. These regulations state that electrical and electronic waste must be separately collected for the proper treatment and recovery to ensure that WEEE is reused or recycled safely. In line with that commitment Natus Medical Incorporated may pass along the obligation for take back and recycling to the end user unless other arrangements have been made. Please contact Natus for details on the collection and recovery systems available to you in your region at natus.com.

Electrical and electronic equipment (EEE) contains materials, components and substances that may be hazardous and present a risk to human health and the environment when WEEE is not handled correctly. Therefore, end users also have a role to play in ensuring that WEEE is reused and recycled safely. Users of electrical and electronic equipment must not discard WEEE together with other wastes. Users must use the municipal collection schemes or the producer/importers take-back obligation or licensed waste carriers to reduce adverse environmental impacts in connection with disposal of waste electrical and electronic equipment and to increase opportunities for reuse, recycling and recovery of waste electrical and electronic equipment.

Equipment marked with the crossed-out wheeled bin is electrical and electronic equipment. The crossed-out wheeled bin symbol indicates that waste electrical and electronic equipment should not be discarded together with unseparated waste but must be collected separately.

Disclaimer:

Natus Medical Incorporated is not responsible for injury, infection or other damage resulting from the use of this product. Any serious incident that has occurred in relation to the device should be reported to Natus Medical Incorporated and the competent authority of the Member State in which the user and/or patient is established. Refer to the Natus website for an electronic copy of this document.

Maintenance

No special maintenance of the XactTrace® belts is required.

Access the eIFU

Instructions for Use are subject to change. The Instructions for Use for this product are available to view and print from Natus.com. View and print the Instructions for Use using Adobe® PDF Reader®. Download a free copy of Adobe PDF Reader from Adobe Systems® at adobe.com. Request a free paper copy of the Instructions for Use in any of these languages by contacting Natus Technical Support or by contacting your local Natus Distributor/Representative.

Glossary of Symbols:

Symbol	Standard Reference	Standard Title of Symbol	Symbol Title as per Referenced Standard	Explanation
	ISO 15223-1 Symbol 5.7.7	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	An indication of Medical Device	This product is a medical device.
	Medical Device Directive 93/42/EEC or Medical Device Regulation 2017/745 as applicable	Council Directive 93/42/EEC Regulation (EU) 2017/745	CE Marking of Conformity	Signifies European technical conformity. Notified body number appears under symbol if applicable.
	ISO 15223-1 Symbol 5.1.2	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Authorized representative in the European Community	Indicates the Authorized representative in the European Community.
	(www.gov.uk/guidance)	UKCA Mark	UKCA Mark	Signifies Great Britain (England, Wales, and Scotland) conformity. UK Approved Body number appears under symbol.
	Swiss Medical Device Ordinance (MedDO)	Swiss Authorized Representative	Swiss authorized representative	Indicates the Authorized representative in Switzerland.
Rx only	21 CFR 801.109(b)(1)	Labeling-Prescription devices	Device is cleared for the US market as requiring a prescription	Indicates that the product is authorized for sale by or on the order of a licensed healthcare practitioner.
	ISO 15223-1 Symbol 5.1.5	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Batch code or Lot code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Symbol 5.1.6	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Symbol 5.1.11	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Country of Origin	Indicates the country of origin.
	ISO 15223-1 Symbol 5.1.1	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Legal Manufacturer	Indicates the medical device manufacturer.
	ISO 15223-1 Symbol 5.1.3	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Symbol 5.3.7	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Temperature limit	Indicates the (storage) temperature limits to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.3.8	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Humidity limit	Indicates the range of (storage) humidity to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.3.9	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Atmospheric pressure limitation	To indicate the acceptable upper and lower limits of atmospheric pressure for transport and storage.
	ISO 15223-1 Symbol 5.4.2	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Do not re-use	Indicates that the medical device is intended for one use, or for use on a single patient during a single procedure.
	WEEE Directive 2012/19/EU	Waste Electrical and Electronic Equipment (WEEE)	Disposal at end of operating life instructions	Indicates that electrical and electronic equipment waste should not be discarded together with unseparated waste but must be collected separately.
	N/A	N/A	Quantity	Number of parts in a package.
	ISO 15223-1 Symbol 5.4.4	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 60601-1 Table D.1 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance		
	Based on ISO 15223-1, Annex B Negation of Symbol 5.4.5	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	Not made with Natural Rubber Latex	Indicates a medical device that is not made with natural rubber latex.
	ISO 60601-1 Table D.2 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance	Follow Instructions for Use	Refer to instruction manual/Booklet. NOTE on ME EQUIPMENT "Follow Instructions for use".
	IEC 60601-1 Table D.2 #2	Medical electrical equipment — Part 1: General requirements for basic safety and Essential performance	Warning	Indicates a hazard of potential personal injury to patient or operator.
	ISO 15223-1 Symbol 5.1.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Importer	Indicates the entity importing the medical device into the locale. This symbol shall be accompanied by the name and address of the importing entity, adjacent to the symbol.
	ISO 15223-1 Symbol 5.1.9	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Distributor	Indicates the entity distributing the medical device into the locale. This symbol shall be accompanied by the name and address of the distributing entity, adjacent to the symbol.
	ISO 15223-1 Symbol 5.1.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied. Use-by date.	Use-by date	Indicates the date after which the medical device is not to be used.

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