Instructions for Use—English

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Indications for Use

Nonin Model 8000SX-Series Reusable Soft Sensor's are indicated for non-invasive pulse-checking and/or continuous monitoring of adult and pediatric patients who are well or poorly perfused. SpO₂ can be used in environments including operating room, surgical recovery, critical care, emergency room, long-term care, home use, and mobile environments, notolament.

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Warnings:

- The use of sensor and oximeter combinations other than Nonin-branded products have not been tested for accuracy as a system and may affect performance of the system. Refer to Nonin pulse oximeter operator's manuals for a complete listing of Nonin-branded oximeters, sensors, and accessories.
- Inspect the sensor application site at least every 6 to 8 hours to ensure correct sensor alignment and skin integrity in the patient. Sensitivity to sensors may vary due to medical status or skin condition.
- Avoid excessive pressure to the sensor application site as this may cause pressure to the skin beneath the sensor.

Cautions:

- Do not use a damaged sensor. If the sensor is damaged, discontinue use immediately.
- Do not sterilize, autoclave or immerse in liquid of any kind.
- Do not use caustic or abrasive cleaning agents on the sensor.
- Follow local governing ordinances and recycling instructions regarding disposal or recycling of the sensor and any components.
- A functional tester cannot be used to assess the accuracy of a pulse oximeter monitor or sensor.
- As with all medical devices, carefully review patient cables and connectors to reduce the possibility of entanglement or strangulation.
- Refer to the pulse oximeter operator's manual for additional warnings and cautions.
- Factors that may degrade pulse oximeter performance include the following:
 - excessive ambient light
 - excessive motion
 - electromagnetic interference
 - moisture in the sensor
 - improperly applied sensor
 - Carboxyhemoglobin
 - methemoglobinemia
 - blood flow restrictors (arterial catheters, blood pressure cuffs, etc.)
 - incorrect sensor type

Symbol	Definition of Symbol
	Follow Instructions for Use
	CAUTION!
	CE Marking indicating conformance to EC Directive No. 93/42/EEC concerning medical devices
	Lot Number
	Protected against vertically falling water droplets (water resistance is tested up to 15 degrees of water spray and solid foreign objects greater than or equal to 2.5 mm in diameter per IEC 60529)
	Storage/Shipping temperature range (if applicable)
	Storage/Shipping humidity range (if applicable)
	Medical prescription required

Choosing the Appropriate Sensor
Use the measurements provided below to determine which sensor should be used. Sensor recommendations are based on digit height (inch/mm), as indicated on the left.

- For heights between 0.5 and 1.0 in. (12.5 – 25.5 mm), use the Model 8000SX (Large).
- For heights between 0.4 and 0.75 inches (10 – 19 mm), use the Model 8000SX (Medium).
- For heights between 0.3 and 0.5 inches (7.5 – 12.5 mm), use the Model 8000SX (Small).

Attaching the Sensor

- Insert the selected digit (refer to the sizing recommendations above) into the sensor as illustrated in Figures 1 and 2. The patient's digit must remain straight.
- Direct the cable along the patient's finger and parallel to the arm leg. (Optional: Secure the sensor cable as needed.)
- Connect the sensor cable to the pulse oximeter or to the patient cable.
- Verify proper operation as described in the pulse oximeter operator's manual.

Note: Proper sensor placement is critical for good measurement. If the sensor is not positioned properly, light may pass the tissue and result in SpO₂ inaccuracies.

Note: The 8000SX-WO2 sensor is compatible with the WristOx2, Model 3150. It is also compatible with Nonin's Model 3100 and 4100 oximeters when used with the 3150W1 adapter.

Cleaning the Sensors

- Cautions:**
 - Clean the sensor before applying to a new patient.
 - Unplug the sensor from the pulse oximeter before cleaning.
 - Do not sterilize, autoclave or immerse the sensor in liquid of any kind. Do not pour or spray any liquids onto the sensor.
 - Do not use caustic or abrasive cleaning agents on the sensor. Do not use cleaning agents containing ammonium chloride. Use of these cleaning agents can shorten the life of the product or void the warranty.
- To clean the sensor, wipe all patient contact surfaces with a soft cloth dampened with a mild detergent or a 10% bleach/90% water solution (household bleach containing less than 10% sodium hypochlorite). Reference sensor in Figure 3.
- Allow the sensor to dry thoroughly before reusing.
- Note:** To minimize cable deterioration when cleaning the cable, gently wipe away the plug end towards the sensor end.

Specifications		
SpO ₂ Accuracy (Adults/Peds) ^{1, 2} :		
Range	Oxygen Saturation (A _{rms} ²) (figure A)	Motion Oxygen Saturation (A _{rms} ²) (figure B)
70 – 100%	±2	±3
70 – 80%	±2	±3
80 – 90%	±2	±2
90 – 100%	±2	±2

SpO₂ Low Perfusion Accuracy: 70% to 100% ±2 digits (A_{max})¹

Pulse Rate Accuracy: 18 to 300 BPM ±3 digits (A_{max})¹

Pulse Rate Low Perfusion Accuracy: 40 to 240 BPM ±3 digits (A_{max})¹

Temperature: 3.4

Operating Temperature: -20 °C to 50 °C (-4 °F to 154 °F)

Humidity: 3.4

Operating Humidity: 10% to 95% non-condensing

Storage/Transportation: 10% to 95% non-condensing

¹A_{max} encompasses 68% of the population.

²Additional accuracy and performance information can be found in the user's manual.

³Accuracy specifications based on Nonin's PureSAT[®] SpO₂ technology and PureSAT[®] PulseLight[®] technology.

⁴For combined oximeter/sensor specifications, refer to the applicable user's manual.

⁵Range as tested with Nonin's PureSAT[®] SpO₂ technology.

Compliance
This product complies with ISO 10993-1.
Not made with natural rubber latex.

Warranty

2 years from the date of delivery.

The device's expected service life is 2 years.

Nonin reserves the right to make changes and improvements to this Instructions for Use and the product it describes at any time, without notice or obligation.

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Nonin Medical, Inc.
13700 1st Avenue North
Plymouth, MN 55441-5443 USA

+1 (763) 553-9968
(800) 356-8874 (US and Canada)
+31 (0)13 - 79 99 040 (Europe)

Fax: +1 (763) 553-7807
+31 (0)13 - 79 99 042 (Europe)

E-mail: info@nonin.com
info@nonin.com (Europe)

Authorized EC Representative:
MPS, Medical Product Service GmbH
Borgasse 20
D-35619 Braunsfeld, Germany

©2016 Nonin Medical, Inc.

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Modelli 8000SX, 8000SX-WO e 8000

Instructions for Use—English

Model 8000SX, 8000SX-WO, 8000SX-WO2 Reusable Soft Pulse Oximeter Sensors

Indications for Use

Nonin's Model 8000SX Series Reusable Soft Sensors are indicated for non-invasive spot-checking and/or continuous monitoring of adult and pediatric patients who are well or poorly perfused. It is intended for use in environments including operating room, surgical recovery, critical care, emergency room, long-term care, home use, and mobile environments.

⚠️ CAUTION: Federal law (USA) restricts this device to sale by or on the order of a licensed practitioner.

Contraindications:

- Do not use the device in an MR environment or in an explosive atmosphere.
- This device is not defibrillation proof per IEC 60601-1.

Warnings:

- The use of sensor and oximeter combinations other than Nonin-branded products have not been tested for accuracy as a system and may affect performance of the system. Refer to Nonin pulse oximeter, operator's manuals for a complete listing of Nonin-branded oximeters, sensors, and accessories.
- Inspect the sensor application site at least every 6 to 8 hours to ensure correct sensor alignment and skin integrity. Patient sensitivity to sensors may vary due to medical status or skin condition.
- Avoid excessive pressure to the sensor application site as this may cause damage to the skin beneath the sensor.

⚠️ Cautions:

- Do not use a damaged sensor. If the sensor is damaged, discontinue use immediately.
 - Do not sterilize, autoclave or immerse in liquid of any kind.
 - Do not use caustic or abrasive cleaning agents on the sensor.
 - Follow local governing ordinances and recycling instructions regarding disposal or recycling of the sensor and any components.
 - A functional tester cannot be used to assess the accuracy of a pulse oximeter monitor or sensor.
 - As with all medical equipment, carefully route patient cables and connections to reduce the possibility of entanglement or strangulation.
 - Refer to the pulse oximeter operator's manual for additional warnings and cautions.
- Factors that may degrade pulse oximeter performance include the following:
- excessive ambient light
 - excessive motion
 - electrostatic interference
 - moisture in the sensor
 - improperly applied sensor
 - Carboxyhemoglobin
 - Methemoglobin
 - blood flow restrictors (arterial catheters, blood pressure cuffs, infusion lines, etc.)
 - incorrect sensor type
 - poor pulse quality
 - venous pulsations
 - anemia or low hemoglobin concentrations
 - carotidocavicular dyes
 - dysfunctional hemoglobin
 - artificial nails or fingernail polish
 - residue (e.g., dried blood, dirt, grease, oil) in the light path

Symbol:	Definition of Symbol
	Follow Instructions for Use
	CAUTION!
	CE Marking indicating conformance to EC Directive No. 93/42/EEC concerning medical devices
	Lot Number
	Protected against vertically falling water drops when enclosure is lifted up to 15 degrees and ingress of solid foreign objects greater than or equal to 2.5 mm in diameter per IEC 60529.
	Storage/shipping temperature range (if applicable)
	Storage/Shipping humidity range (if applicable)
	Medical prescription required

Choosing the Appropriate Sensor

- Use the measurements provided below to determine which sensor should be used. Sensor recommendations are based on digit height (thickness), as indicated at left.
- For heights between 0.5 and 1.0 in. (12.5 – 25.5 mm), use the Model 8000SL (Large).
- For heights between 0.4 and 0.75 inches (10 – 19 mm), use the Model 8000SM (Medium).
- For heights between 0.3 and 0.5 inches (7.5 – 12.5 mm), use the Model 8000SS (Small).

Attaching the Sensor

- Insert the selected digit (refer to the sizing recommendations above) into the sensor as illustrated in Figures 1 and 2. The patient's digit must reach the end of the sensor.
- Direct the cable along the patient's finger/toe, parallel to the arm/leg. *(Optional: Secure the sensor cable as needed.)*
- Connect the sensor cable to the pulse oximeter or to the patient cable.
- Verify proper operation as described in the pulse oximeter operator's manual.

Note: Proper sensor placement is critical for good performance. If the sensor is not positioned properly, light may bypass the tissue and result in SpO₂ inaccuracies.

Note: The 8000SX-WO2 sensor is compatible with the WristOx₂ Model 3150. It is also compatible with Nonin's Model 3100 and 4100 oximeters when used with the 3150WI adapter.

Cleaning the Sensors

⚠️ Cautions:

- Clean the sensor before applying it to a new patient.
- Unplug the sensor from the pulse oximeter before cleaning.
- Do not sterilize, autoclave or immerse the sensor in liquid of any kind. Do not pour or spray any liquids onto the sensor.
- Do not use caustic or abrasive cleaning agents on the sensor. Do not use cleaning agents containing ammonium chloride. Use of these chemicals may shorten the life of the product.

- To clean the sensor, wipe all patient contact surfaces with a soft cloth dampened with a mild detergent or a 10% bleach/90% water solution (household bleach [containing less than 10% sodium hypochlorite]). *Reference sensor in Figure 3.*
- Allow the sensor to dry thoroughly before reusing.

Note: To minimize cable deterioration when cleaning the cable, gently wipe away from the plug end towards the sensor end.

Specifications	
SpO ₂ Accuracy (Adults/Peds) ^{1, 2} :	
Range	Motion Oxygen Saturation (A _{mes} *) (figure A)
70 – 100%	±2
70 – 80%	±3
80 – 90%	±2
90 – 100%	±2
SpO ₂ Low Perfusion Accuracy:	
70% to 100% ±2 digits (A _{mes}) ¹	
Pulse Rate Accuracy:	
18 to 300 BPM ±3 digits (A _{mes}) ¹	
Pulse Rate Low Perfusion Accuracy:	
40 to 240 BPM ±3 digits (A _{mes}) ¹	
Temperature:	
-20 °C to 50 °C (-4 °F to 122 °F)	
Storage/Transportation:	
-40 °C to 70 °C (-40 °F to 158 °F)	
Humidity:	
10% to 95% non-condensing	
Storage/Transportation:	
10% to 95% non-condensing	
* ±1 A _{mes} encompasses 68% of the population.	
† Additional accuracy and performance information can be found in the sensor accuracy document on the operator's manual CD.	
‡ Accuracy specifications based on Nonin's PureSAT™ SpO ₂ technology and Puret-gH [®] sensor technology.	
§ For combined oximeter/sensor specifications, refer to the applicable oximetry system's operator's manual.	
¶ Gasas as tested with Nonin's PureSAT SpO ₂ technology.	
Measurement Wavelengths and Output Power**	
Red:	660 nanometers @ 0.8 mW nominal
Infrared:	910 nanometers @ 1.2 mW nominal
** This information is especially useful for clinicians performing photodynamic therapy.	
Compliance	
This product complies with ISO 10993-1.	
Not made with natural rubber latex.	
Warranty	
2 years from the date of delivery.	
The device's expected service life is 2 years.	
Nonin reserves the right to make changes and improvements to this Instructions for Use and the product it describes at anytime, without notice or obligation.	

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Indications for Use

Nonin's Model 8000SX Series Reusable Soft Pulse Oximeters are indicated for non-invasive spot-checking and/or continuous monitoring of adult and pediatric patients who are well or poorly perfused. It is intended for use in environments including operating room, pre/post recovery, critical care, emergency room, long-term care, home use, and mobile environments.

Key CAUTION: Federal law (USA) restricts this device to sale by or on the order of a licensed practitioner.

Contraindications:

- Do not use the device in an MR environment or in an explosive atmosphere.
- This device is not defibrillation proof (IEC 60601-1).

Warnings:

- The use of sensor and oximeter combinations other than Nonin-branded products have not been tested for accuracy as a system and may affect performance of the system. Refer to Nonin pulse oximeter operator's manuals for a complete listing of Nonin-branded oximeters, sensors, and accessories.
- Inspect the sensor application site at least every 8 to 8 hours to ensure correct sensor alignment and skin integrity. Patient sensitivity to sensors may vary due to medical status or skin condition.
- Avoid excessive pressure to the sensor application site as this may cause damage to the skin beneath the sensor.

Cautions:

- Do not use a damaged sensor. If the sensor is damaged, discontinue use immediately.
- Do not sterilize, autoclave or immerse in liquid of any kind.
- Do not use caustic or abrasive cleaning agents on the sensor.
- Follow local governing ordinances and recycling instructions regarding disposal or recycling of the sensor and any components.
- A functional tester cannot be used to assess the accuracy of a pulse oximeter monitor or sensor.
- As with all medical equipment, carefully route patient cables and connections to reduce the possibility of entanglement or strangulation.
- Refer to the pulse oximeter operator's manual for additional warnings and cautions.
- Factors that may degrade pulse oximeter performance include the following:
 - excessive ambient light
 - excessive motion
 - electromagnetic interference
 - moisture in the sensor
 - prolonged applied sensor
 - Carboxyhemoglobin
 - Methemoglobin
 - blood flow restrictors (arterial catheters, blood pressure cuffs, restraint, etc.) in the right path
 - residue (e.g. dried blood, dirt, infusion lines, etc.)
 - incorrect sensor type
 - poor pulse quality
 - venous pulsations
 - anemia or low hemoglobin concentrations
 - disfunctional hemoglobin
 - artificial nails or fingernail polish
 - residue (e.g. dried blood, dirt, infusion lines, etc.)

Symbols:

Symbol	Definition of Symbol
	Follow Instructions for Use
	CAUTION!
	CE Marking indicating conformance to EC Directive No. 93/42/EEC concerning medical devices
	Lot Number
	Protected against vertically falling water drops when enclosure is tilted up to 15 degrees and ingress of solid foreign objects greater than or equal to 2.5 mm in diameter per IEC 60529.
	Storage/shipping temperature range (if applicable)
	Storage/shipping humidity range (if applicable)
	Medical prescription required

Choosing the Appropriate Sensor

- Use the following dimensions provided below to determine which sensor should be used. Sensor recommendations are based on digit height (thickness), as indicated at left.
- For heights between 0.5 and 1.0 in. (12.5 – 25.5 mm), use the Model 8000SL (Large).
- For heights between 0.4 and 0.75 inches (10 – 19 mm), use the Model 8000SM (Medium).
- For heights between 0.3 and 0.5 inches (7.5 – 12.5 mm), use the Model 8000SS (Small).

Attaching the Sensor

- Insert the selected digit (refer to the sizing recommendations above) into the sensor as illustrated in Figures 1 and 2. The patient's digit must reach the end of the sensor.
- Direct the cable along the patient's finger, parallel to the artery. (Optional: Secure the sensor cable as directed below.)
- Connect the sensor cable to the pulse oximeter or to the patient cable.
- Verify proper operation as described in the pulse oximeter operator's manual.
- Note: Proper sensor placement is critical for good performance. If the sensor is not positioned properly, light may bypass the tissue and result in SpO₂ inaccuracies.
- Note: The 8000SX-WO2 sensor is compatible with the WristOx[®] 3150. It is also compatible with Nonin's Model 3100 and 4100 oximeters when used with the 3150W1 adaptor.

Cleaning the Sensors

- Cautions:**
 - Clean the sensor before applying it to a new patient.
 - Unplug the sensor from the pulse oximeter before cleaning.
 - Do not sterilize, autoclave or immerse the sensor in liquid of any kind. Do not pour or spray any liquids onto the sensor.
 - Do not use caustic or abrasive cleaning agents on the sensor. Do not use cleaning agents containing ammonium chloride. Use of these chemicals may shorten the life of the product.
- Clean the sensor using a soft cloth dampened with a mild detergent or a 10% bleach/90% water solution (household bleach containing less than 10% sodium hypochlorite). Reference sensor in Figure 3.
- Note: To minimize cable deterioration when cleaning the cable, gently wipe away from the plug end towards the sensor end.

Specifications

Range	Oxygen Saturation (A _{max}) ¹ (figure A)	Motion Oxygen Saturation (A _{max}) ¹ (figure B)
70 – 100%	±2	±3
70 – 80%	±2	±3
80 – 90%	±2	±2
90 – 100%	±2	±2

SpO₂ Low Perfusion Accuracy: 70% to 100% ±2 digits (A_{max})¹
Pulse Rate Accuracy: 18 to 300 BPM ±3 digits (A_{max})¹
Pulse Rate Low Perfusion Accuracy: 40 to 240 BPM ±3 digits (A_{max})¹

Temperature: 3.4
Operating: -20 °C to 50 °C (-4 °F to 122 °F)
Storage/Transportation: -40 °C to 70 °C (-40 °F to 158 °F)
Humidity: 3.4
Operating: 10% to 95% non-condensing
Storage/Transportation: 10% to 95% non-condensing

***1** A_{max} encompasses 68% of the population.
***2** Accuracy specifications based on Nonin's PureSat[®] SpO₂ technology and PureLight[®] sensor technology.
***3** For combined oximeter/sensor specifications, refer to the applicable oximetry system's operator's manual.
***4** As tested with Nonin's PureSAT SpO₂ technology.

Measurement Wavelengths and Output Power**

Red: 660 nanometers @ 0.8 mW nominal
Infrared: 910 nanometers @ 1.2 mW nominal
** This information is especially useful for clinicians performing photodynamic therapy.

Compliance

This product complies with ISO 10993-1.
Not made with natural rubber latex.

Warranty

Nonin warrants the device from the date of delivery.
The device's expected service life is 2 years.
Nonin reserves the right to make changes and improvements to these Instructions for Use and the product it describes at any time, without notice or obligation.

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