
LAST 3 REVS	REV	DESCRIPTION	DCO/ECO	BY	APPR	DATE
	01A	Update drawing to include correct size formatting for printing	ECO#21093	MPR	MES	06-APR-16
	01B	Update to correct typo - "rericts" should read "restricts"	ECO#21946	MPR	MES	22-JUN-16
	02	Update INTENDED USE STATEMENT	ECO#22498	MES	JT	20-SEP-16

Y RoHS Compliant to 2011/65/EU Y Latex Free Y BPA Free Y DEHP <0.1% by Weight


Leaflet size 830 mm x 210 mm (flat)
10 page folded (page size 83 mm)

See pages 2 & 3 for IFU content

NOTES.

- 1. INSTRUCTIONS FOR USE; LANGUAGE TRANSLATIONS REQUIRED:- ENGLISH, FRENCH, GERMAN, SPANISH, FINNISH, DANISH, SWEDISH, NORWEGIAN, POLISH, DUTCH, GREEK, ITALIAN, CZECH, HUNGARIAN, TURKISH, PORTUGESE, CHINESE, JAPANESE, RUSSIAN.
- 2. LEAFLET TO BE DOUBLE SIDED PRINT. NINE LANGUAGES PAGE 1, TEN LANGUAGES PAGE 2, 19 languages.
- 3. PRINT COLOUR: ALL TEXT & GRAPHICS, BLACK.
- 4. OUTLINE IS FOR REFERENCE ONLY; DO NOT PRINT IT.
- 5. FINISHED LEAFLET TO BE SUPPLIED FAN FOLDED, EACH PAGE WIDTH AS STATED ABOVE. COMPANY DETAILS TO FACE OUTWARDS.
- 6. PRINTING PAPER TO BE BRIGHT WHITE, APPROXIMATELY 80 GSM.


THIS DOCUMENT CONTAINS CONFIDENTIAL OR PROPRIETARY INFORMATION OF NATUS NEUROLOGY INCORPORATED. NEITHER THIS DOCUMENT NOR THE INFORMATION HEREIN IS TO BE REPRODUCED, DISTRIBUTED, USED OR DISCLOSED, IN PART OR IN WHOLE, EXCEPT AS SPECIFICALLY AUTHORIZED BY N	UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES ARE: .X± .XX± .XXX± REMOVE ALL BURRS, BREAK SHARP EDGES DRAWING PER ANSI Y14.5 WELDING SYMBOLS PER AWS A2.0 DO NOT SCALE DRAWING		MATERIAL See notes		natus neurology Natus Manufacturing Limited IDA Business Park Gort, Co.Galway, Ireland		
			FINISH See notes				
	DESIGN BY Mary Ryan		DATE 01-Mar-2016		TITLE IFU for Boject DHN		
	APPROVED BY Mike Scanlon		DATE 01-Mar-2016		SIZE B	DRAWING/PART NO. 017669	REV 02
					SCALE: NTS E-FILE:		SHEET 1 OF 3

Bo-ject® DHN
Disposable Hypodermic
Needle Electrode

Instructions for Use

English
Intended use
The disposable hypodermic monopolar needles are FOR SINGLE USE ONLY in electromyography (EMG) in situations wherein it is desired to insert an electrode, in the form of a probe, into a patient to locate a particular muscle and then inject a medication into that muscle.

Instructions for use
Remove the BoJect® DHN electrode from its outer package.
Syringes with Luer lock: Push syringe firmly onto electrode and turn clockwise until devices lock. Luer slip syringes may be pushed firmly into electrode and secured with a slight downward rotation movement.


Connection Instructions for all other Electrode Systems:
Plug the needle into the active recording or stimulating connections of your instrument. Connect the necessary reference and ground electrodes as appropriate for stimulating or recording.
Check expiration date on package before use.
Usage: Needle Electrodes are for professional use only.

Caution:
Federal law (USA), restricts this device to sale by or on the order of a licensed medical practitioner.

Table with 2 columns: Symbol, English. Rows include Needle length and gauge, Latex status, Phthalates status, Federal law (USA) restrictions, Single use only, Instructions for use, Do not use if package is damaged, Caution/Warning, and Sterility.

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
Table with 2 columns: Symbol, Czech. Rows include Needle length and gauge, Latex status, Phthalates status, Federal law (USA) restrictions, Single use only, Instructions for use, Do not use if package is damaged, Caution/Warning, and Sterility.

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
Table with 2 columns: Symbol, Dansk. Rows include Needle length and gauge, Latex status, Phthalates status, Federal law (USA) restrictions, Single use only, Instructions for use, Do not use if package is damaged, Caution/Warning, and Sterility.

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
Table with 2 columns: Symbol, German. Rows include Needle length and gauge, Latex status, Phthalates status, Federal law (USA) restrictions, Single use only, Instructions for use, Do not use if package is damaged, Caution/Warning, and Sterility.

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
Table with 2 columns: Symbol, Greek. Rows include Needle length and gauge, Latex status, Phthalates status, Federal law (USA) restrictions, Single use only, Instructions for use, Do not use if package is damaged, Caution/Warning, and Sterility.

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
Table with 2 columns: Symbol, Finnish. Rows include Needle length and gauge, Latex status, Phthalates status, Federal law (USA) restrictions, Single use only, Instructions for use, Do not use if package is damaged, Caution/Warning, and Sterility.

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
Table with 2 columns: Symbol, French. Rows include Needle length and gauge, Latex status, Phthalates status, Federal law (USA) restrictions, Single use only, Instructions for use, Do not use if package is damaged, Caution/Warning, and Sterility.

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
Table with 2 columns: Symbol, Magyar. Rows include Needle length and gauge, Latex status, Phthalates status, Federal law (USA) restrictions, Single use only, Instructions for use, Do not use if package is damaged, Caution/Warning, and Sterility.

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Table with 2 columns: Symbol, Italian. Rows include Needle length and gauge, Latex status, Phthalates status, Federal law (USA) restrictions, Single use only, Instructions for use, Do not use if package is damaged, Caution/Warning, and Sterility.

Natus Manufacturing Limited
IDA Business Park
Gort, Co. Galway, Ireland
Tel: +1-800-356 0007 or +1-608-829-8500
Fax: +1-608-829-8516
supplies.nc@natus.com

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