REGULATORY: Class I device. Complies with MDR (EU) 2017/745 and manufactured according to the procedures of BIONEN Quality System certified ISO 13485.

GUARANTEE AND LIMITATION OF LIABILITY:

BIONEN guarantees that the device is free of material and manufacturing defects for a period of twelve months from the date of purchase. The sole responsibility of BIONEN and our dealers is limited to the replacement or repair of the product, in the opinion of BIONEN, without charge for parts or work, if it is proven that one of the parts is defective in manufacturing, performance or materials, during warranty period. Under no circumstances will BIONEN or its dealers be liable for any loss of earnings or for damages, direct, indirect or accidental, including loss of profit, damage to property, or personal injury arising from the use, or from the inability to use the product. This warranty only covers the original purchaser and supersedes all other warranties or prior agreements, either expressed or implied. This warranty is void if the product is used for purposes other than those for which it is intended, or if it is subject to abuse, misuse, tampering, negligence, or unauthorized changes. Use of this product constitutes acceptance of the full warranty. Notify the BIONEN Quality Assurance Service for any malfunction or defect found on this product.

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MUSCLE STIMULATOR

REF.	VERSIONS
0038.0xx	Flat bridge stimulator with belt / handle
0038.02x 0038.03x	Flat bridge stimulator with belt and pads
0038.2xx	Stimulator with velcro fastening
0039.0xx	Bipolar muscle stimulator
0039.20x	Monopolar muscle stimulator

INTENDED USE: EMG - EP

INSTRUCTION FOR USE: Lean the stimulator's electrode on the point to be stimulated, applying light pressure.

In case of stimulators with pads, dip the pads in a 0,9% physiological saline solution. Clean up patient's skin, then apply the stimulator using the belt.

Connect the plugs to the device, making sure to reach the best contact impedance.

CLEANING INSTRUCTIONS: Carefully clean the stimulator after each use with a damp cloth or gauze. Do not immerse in any liquid.

For stimulators with pads, clean the felts after each test with water and a neutral detergent.

DISINFECTION: Disinfect the device using either Amocid, Lysoformin, Alhyndex, IndicinGG or equivalent products, even at high concentrations, without immersing the device into the liquids.

STERILIZATION METHOD: NON-STERILE device that does not require sterilization before use.

LIFETIME: After 5 years from sale, is not guaranteed the proper functionality of the device.

STORAGE:

No special storage precautions are required. However, it is recommended to store the electrodes in their original box after cleaning.

CONTRINDICATIONS AND INTERACTIONS: none

POSSIBLE TOXICITY DECLARED: none. LATEX FREE

WARNINGS AND PRECAUTIONS: For professional use only.

TRANSPORT AND DISPOSAL

No special transportation or disposal requirements.

HANDLING AND STORAGE: Do not overlap excessive loads.

WASTE DISPOSAL: Manage the product as sanitary waste in conformity with the regulations in force.

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