Disposable Subdermal Needle Electrode

Instructions for use

Agoelettrodo subdermale monouso

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Spes Medica S.r.l.
Via Buccari 21
16153 Genova (GE), Italy
Tel: +39 010 390343
Fax: +39 010 3072345
Web: www.spesmedica.com

ENGLISH

Disposable Subdermal Needle Electrode

Intended use

Spes Medica Subdermal Needle Electrodes are intended for use with recording, monitoring and stimulation equipment for the purpose of recording of biopotential signals. Examples include: Electromyography (EMG), Electroencephalography (EEG) and Nerve potential signals. The electrodes are sterile and for single patient use only.

Instructions for use

- Do not use the selected needle if the pouch is damaged;
- Open the pouch and remove needle from sterile pouch;
- Connect the needle with the equipment;
- Remove the needle sheath;
- Introduce the needle into the area to be tested. Insert the needle up to hub.
- After the use, dispose the needle according to standard protocols or alternatively put the needle in approve biohazard sharps containers.

The electrodes are sterile and for single patient only

Warnings

- The metal part of the needle electrode mustn't be visible after insert. Make sure of the total insert of the needle electrode into the patient.
- The subdermal needle electrode is for single patient only.
- This single-use product is not designed for re-use. Attempts to clean and re-use this single-use device expose patients and operators to a risk of cross contamination.
- Do not re-sterilize
- Not compatible with MR
- In case of subdermal needles or cork screws be used together with high voltage and frequency instruments (e.g. electrosurgery), a minimum distance of 20-30 cm between subdermal needles cables and stimulation instrument cables is required otherwise there could be electrical interferences leading to skin burns.
- Maximum stimulation value: 60 mA.

	Complies with the EC directive 93/42/EEC			
	Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices			
	Catalogue number		Lot number	
	Single Use Only		Do not sterilize	
	Consult instruction for use		Do not use if package is damaged	
	Warnings and Cautions associated with this device.			Expire by
	Ethylene Oxide Sterilization			MR Unsafe
	Manufacturer			

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