Dr. Langer Medical GmbH

Am Bruckwald 26
79183 Waldkirch
GERMANY

Fon: +49 (0) 76 81 / 47 45 4-0
Fax: +49 (0) 76 81 / 47 45 4-29

www.medical-langer.de
info@medical-langer.de

DR. LANGER
M E D I C A L

EINMAL-STIMULATIONSSONDEN, Klasse III Gebrauchsanweisung	DEUTSCH	2
DISPOSABLE STIMULATION PROBES, Class III Instructions for Use	ENGLISH	8
SONDES DE STIMULATION À USAGE UNIQUE, classe III Notice d'utilisation	FRANÇAIS	14
SONDE DI STIMOLAZIONE MONOUSO, classe III Istruzioni per l'uso	ITALIANO	20
SONDAS DE ESTIMULACIÓN DESECHABLES, clase III Instrucciones de uso	ESPAÑOL	26
SONDAS DE ESTIMULAÇÃO DE USO ÚNICO, classe III Instruções de utilização	PORTUGUÊS	32
WEGWERPSTIMULATIESSONDEN, Klasse III Gebruiksaanwijzing	NEDERLANDS	38
JEDNORAZOWE SONDY DO STYMULACJI, Klasa III Instrukcja obsługi	POLSKI	44
JEDNORÁZOVÉ STIMULAČNÍ ELEKTRODY, třída III Návod k použití	ČESKY	50
ΝΕΥΡΟΔΙΕΓΕΡΤΕΣ ΜΙΑΣ ΧΡΗΣΗΣ, Κατηγορία III Οδηγίες χρήσης	ΕΛΛΗΝΙΚΑ	56
ОДНОРАЗОВЫЕ СТИМУЛИРУЮЩИЕ ЗОНДЫ, Класс III Инструкции по эксплуатации	РУССКИЙ	63
TEK KULLANIMLIK STİMÜLASYON PROBLARI, Sınıf III Kullanım Kılavuzu	TÜRKÇE	70

Dr. Langer Medical GmbH
Am Bruckwald 26
79183 Waldkirch
GERMANY

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Instructions for use Art. No.: 05-0070

DISPOSABLE STIMULATION PROBES

Class III in accordance with Directive 93/42/EEC

Monopolar disposable stimulation probe (Art. No. 40-1303), bipolar disposable microfork probe (Art. No. 40-1304), monopolar disposable ball probe (Art. No. 40-1306), bipolar concentric disposable stimulation probe (Art. No. 40-1307), bipolar disposable microfork ball probe (Art. No. 40-1308)

CONTENTS, ENGLISH

1	Intended use.....	8
2	Safety notices and precautions.....	9
3	Symbols used on the packaging.....	11
4	Combination and connection of the disposable stimulation probes.....	11
5	Handling	13
6	Troubleshooting.....	13
7	Disposal.....	13
8	Technical specifications.....	13
9	Warranty	13
10	Copyright and reproduction of the instructions for use.....	13

1 Intended use

1.1 Intended purpose

The disposable stimulation probes of Dr. Langer Medical GmbH are intended for single use in the field of intraoperative neurological monitoring. The devices are disposable stimulation probes for temporary direct stimulation of nerves, muscles and the central nervous system (CNS).

The disposable stimulation probes of Dr. Langer Medical GmbH are used to localise the structures mentioned above and to verify their excitability and integrity. With the help of evoked biopotentials it is possible to confirm the functioning of nerves and muscles, map their pathways or identify specific areas in the brain. Direct electrical stimulation of nerves, muscles or the CNS is performed via the electrodes at the tip of the disposable stimulation probe. The design of the monopolar and bipolar electrode tips prevents tissue injuries when the probes are employed in compliance with their intended use.

Dr. Langer Medical GmbH recommends using their disposable stimulation probes with AVALANCHE® neuromonitors and TWISTER® neurostimulators. All information about connections and technical specifications of neuromonitors and neurostimulators refer to the medical devices AVALANCHE® and TWISTER® of Dr. Langer Medical GmbH.

Dr. Langer Medical GmbH offers the following types of disposable stimulation probes:

Model	Article No.	Illustration
Monopolar probe with insulated shaft	40-1303	
Microfork probe, bipolar	40-1304	
Ball probe, monopolar	40-1306	
Bipolar concentric probe, angled tip	40-1307	
Microfork ball probe, bipolar	40-1308	

1.2 Indication

The Dr. Langer Medical GmbH disposable stimulation probes are used intraoperatively in the following areas of application:

Articles 40-1303, 40-1304, 40-1307, 40-1308

- Direct stimulation of motor nerves, muscles and the central nervous system (CNS) in order to localise them and monitor their integrity during surgical interventions.
- Localisation and monitoring of the integrity of the motor sphincter muscles in the pelvic floor.

Article 40-1306

- Transmission of electrical stimulation pulses to pedicle screws, drilling instruments or drilled holes in order to check the position of the screws in relation to the motor nerves in the spine.
- Direct cortical stimulation and localisation of the motor areas in the brain and the speech centre.

1.3 Contraindication

The disposable stimulation probes are not intended for direct contact with the heart or the central circulatory system. The disposable stimulation probes may not be used in patients with active implants, such as cardiac pacemakers or ICDs, or in pregnant patients.

1.4 User profile

The disposable stimulation probes of Dr. Langer Medical GmbH may only be used in patients for the indication listed in section 1.2 by physicians who meet the following requirements:

- They have been trained in the proper handling of the medical device.
- They have read and understood the instructions for use.
- They are experienced in handling sterile disposable instruments.
- They have an understanding of neurophysiology especially regarding intraoperative nerve and muscle stimulation, stimulation of the CNS and recording of MUAP.

- They are familiar with the areas of application of intraoperative localisation of nerves, muscles and specific areas of the brain and monitoring their integrity, as well as with the applicable contraindications.
- They know how to recognise and prevent false negative stimulation reactions during direct nerve stimulation.
- They know how to safely use electric stimulators.
- They are aware of the effects of muscle relaxants on EMG signals.

The disposable stimulation probes of Dr. Langer Medical GmbH may only be prepared and passed on to the surgeon, at the surgeon's request, by OR staff who meet the following requirements:

- They have been trained in the proper handling of the medical device.
- They have read and understood the instructions for use.
- They are experienced in handling sterile disposable instruments, their application, storage, verification of expiry date and disposal.

In conjunction with neuromonitors and neurostimulators, the disposable stimulation probes of Dr. Langer Medical GmbH are intended to assist the surgeon in locating motor nerves, muscles and the CNS and verifying their integrity during surgical interventions. This, however, is no substitute for experience in surgical interventions and making appropriate decisions.

1.5 Patient profile

In conjunction with neuromonitors and neurostimulators, the disposable stimulation probes of Dr. Langer Medical GmbH are applied during surgical interventions involving the motor nerves and muscles according to the indication in section 1.2. There are no known clinical limitations regarding age or gender of the patients.

1.6 Exclusions from the intended purpose

The administration of muscle relaxants will partially or completely block the stimulus conduction from the motor nerve to the innervated muscle. In this case, the device cannot fulfil its intended purpose, i.e., the localisation and functional check of motor nerves and muscles as well as cortical mapping by means of electrical stimulation.

1.7 Use, transport, storage and disposal

The disposable stimulation probes are intended for intraoperative application in the OR. For information about handling the probes, please refer to chapter 5 in these instructions for use. For detailed information about transport and storage of the disposable stimulation probes, please refer to the device labelling and to chapter 3 in these instructions for use. For information about disposal of the medical device, please refer to chapter 7. The disposable stimulation probes may not be disposed of as unsorted municipal waste. They must be disposed of in compliance with the regulations for the disposal of surgical waste applicable in your facility.

2 Safety notices and precautions

The requisite personal protective equipment must be provided by the operator. All existing safety devices must be checked regularly.

2.1 The meaning of safety notices

 DANGER	• indicates an imminent hazard that can result in death if it is not avoided. • indicates an imminent hazard that can result in devastation of the environment if it is not avoided.
 WARNING	• indicates a potential hazard that can result in serious bodily injuries if it is not avoided. • indicates a potential hazard that can affect the proper functioning of nearby equipment if it is not avoided. • indicates a potential hazard that can result in significant damage of the environment if it is not avoided.
 CAUTION	• indicates a potential hazard that can result in mild or moderate bodily injuries if it is not avoided. • indicates a potential hazard that can result in brief malfunctions of the device if it is not avoided.
ATTENTION	• indicates a potential hazard that can result in property damage if it is not avoided.

2.2 The meaning of notices

"Notice:"

- indicates manufacturer information which directly or indirectly concerns the safety of individuals or the protection of property. The information does not directly concern a hazard or a hazardous situation.
- indicates manufacturer information that is important or useful for the operation or maintenance of the product.

2.3 Knowledge of the instructions for use

Understanding the instructions for use is absolutely essential for correct operation of the product.

Therefore, everyone who

- prepares
- works with
- and cleans/reprocesses

the device, must read the instructions for use. Observe the safety notices of each chapter.

2.4 Informal safety measures

These instructions for use must be readily accessible at the location where the devices are used. The safety notices and precautions must be kept in a legible condition.

2.5 Observing the safety notices

Working with a medical device is generally associated with residual risks to patients, to the medical staff and to the environment. Risks cannot be entirely eliminated by design measures alone. Safety does not depend solely on the product. Safety depends largely on the training of medical staff and correct operation of the product. Everyone working with the device is expected to read, understand and observe the safety notices in this chapter.

2.6 Safety notices and precautions

The safety notices are divided into the following hazard categories: Defective device / defective accessories, application and clinical residual risks, anaesthesia.

2.6.1 Defective device, defective accessories

 CAUTION	Defective/modified device or connected accessories When using a defective/modified device and/or accessories, safe operation and handling cannot be guaranteed. The user and patient may suffer injuries. • Check the device and all accessories for damage before each use (mechanical damage, kinked cable, damaged insulation, bent or broken probe tips...). • A damaged device and/or accessory may not be used. Discard defective devices and accessories.
--	---

2.6.2 Application and clinical residual risks

 WARNING	Use of device from damaged pack/past its expiry date There is a risk of infection when a device from a damaged pack and/or after its expiry date is used. • Check the packaging for damage. • Check the expiry date.
 WARNING	Use of a reprocessed device The disposable stimulation probes are devices intended for single use. They are not intended to be reprocessed. The procedures employed for the reprocessing of the device with the purpose of reusing it, including its repeated sterilisation, may degrade its mechanical and biological characteristics to an unacceptable level and considerably increase the risk of contamination. There is a risk of infection, allergies, toxic reactions and other intolerances. • After use, dispose of the device according to chapter 7 of these instructions for use. • Do not reprocess the device.
 WARNING	Wrong position of reference electrode¹ during monopolar stimulation with respect to the HF neutral electrode The current from the electrosurgical unit is conducted via the reference electrode. Risk of skin burns at the reference electrode. • When performing monopolar stimulation, do not attach the reference electrode in the path of the current of the electrosurgical unit from the surgical site to the HF neutral electrode.
 WARNING	Simultaneous use of an electrosurgical unit and the device, or device in contact with HF electrode. Risk of skin burns at reference electrode, lead electrode or stimulation electrode or at other conductive parts. • Do not allow the device to come into contact with the HF electrode. • Do not activate stimulators for direct nerve / muscle stimulation at the same time as the electrosurgical unit.
 WARNING	Use with neuromonitors/neurostimulators not designed for this purpose Using the devices with neuromonitors/neurostimulators not indicated for the intended purpose of the stimulation probes may cause high current densities at the probe tips. There is a potential risk of causing necrosis. • Use the devices only with suitable neuromonitors/neurostimulators. • The devices must only be employed in compliance with their intended use, see chapter 1.
 WARNING	Muscle does not respond although nerve/muscle is intact (false negative response to stimulation) The stimulation current actually applied to the nerve may be reduced by current division according to Ohm's law, particularly due to fluid accumulated in the area of the electrode-tissue contact area or when the stimulation electrode is too far away from the excitable structure. • Prevent fluid accumulation.
 WARNING	Stimulation not possible Dry tissue and/or dried tissue residue at the probe tip causes a very high contact impedance between probe tip and tissue. It is not possible to deliver a stimulation pulse. • Clean the probe tip. • Moistened the tissue.
 WARNING	Use of neuromonitors/stimulators not intended for direct nerve / muscle stimulation When using the device in combination with electrical neuromonitors/stimulators not intended for direct nerve / muscle stimulation, there is a risk of high current densities at the probe tips. There is a risk of skin burns and tissue lesions. • Use the device only in compliance with its intended use described in chapter 1. • Observe the instructions for use of the neuromonitor/stimulator.
 WARNING	Epileptic seizure induced by direct cortical stimulation Direct cortical stimulation may induce seizures. • Review the selected stimulation parameters and contact the anaesthetist.
 CAUTION	Perforation of blood vessels When using the device for surgical preparation and with excessive pressure, there is a risk of perforating blood vessels. This may cause bleeding and haematomas. • Avoid excessive pressure on the tissue during stimulation. • Review the anatomical structures in the area to be stimulated.
 CAUTION	Needle-stick injuries or infections due to needle electrodes¹ when using monopolar stimulation probes During and after use, there is a risk of needle-stick injuries. Used needles are potentially contaminated. This may cause infections, bleeding and haematomas. • Handle needles with great care. Unnecessary needle-stick injuries must be avoided. • Review the anatomical structures in the area where the needle will be inserted.

¹ When using monopolar stimulation instruments, an extra needle electrode as reference electrode is required. Position the reference electrode in the tissue near the stimulation site. See chapter 4.2.

2.6.3 Anaesthesia

WARNING

Use of muscle relaxants

The administration of muscle relaxants will partially or completely block the stimulus conduction from the motor nerve to the innervated muscle. In this case, the device cannot fulfil its intended purpose, i.e., the localisation and functional control of nerves and muscles by means of electrical stimulation.

- Do not use the device if muscle relaxants are administered.
- Contact the anaesthetist.

CAUTION

Use in non-anaesthetised patients and influence of anaesthetics

Application of the device may cause sensations of pain and lead to muscular responses (twitch).

- Use the device only in anaesthetised patients.
- Contact the anaesthetist.

3 Symbols used on the packaging

	Identifies the manufacturer of the medical device according to the EU Directive 93/42/EEC.		Identifies a single sterile barrier system
	Article number. Indicates the manufacturer's order number which allows the medical device to be identified.		Identifies a single sterile barrier system with protective packaging outside.
	Lot code. Indicates the manufacturer's lot code, which allows the batch or lot to be identified.		Identifies a disposable or single-patient use/one-time-use medical device.
	Indicates the date of manufacture.		Identifies a medical device which does not contain any natural rubber latex.
	Indicates the expiry date after which the use of the medical device is no longer permitted.		Identifies a medical device.
	Indicates a medical device which was sterilised with ethylene oxide.		Unique Device Identification
	Attention. Indicates that the user is required to study the instructions for use, particularly with respect to safety information such as warnings and precautions that, for various reasons, cannot be placed directly on the medical device.		GS1 Datamatrix 2D code 2D code as legible text: (01)04251807 200186 (11)210100 (17)230100 (10)123456
	Observe instructions for use. Indicates that the user is required to refer to the instructions for use.		Symbol indicates compliance with the applicable EU Directives. 0297: ID number of the notified body (DQS Medizinprodukte GmbH, Frankfurt am Main, Germany)
	Indicates a medical device which should not be used if its packaging is damaged or open.		
	Identifies a medical device which must be protected from light sources.		Identifies a medical device which must be protected from humidity.

4 Combination and connection of the disposable stimulation probes

4.1 Combination with accessories and other medical devices

Accessories and combinable medical devices manufactured by Dr. Langer Medical GmbH are recommended for use with the disposable stimulation probes of Dr. Langer Medical GmbH. All information about connections and technical specifications of accessories and medical devices provided in this document refer to products from Dr. Langer Medical GmbH.

WARNING

Use with neuromonitors/neurostimulators not designed for this purpose

Using the devices with neuromonitors/neurostimulators not indicated for the intended purpose of the stimulation probes may cause high current densities at the probe tips. There is a potential risk of causing necrosis.

- Use the devices only with suitable neuromonitors/neurostimulators.
- The devices must only be employed in compliance with their intended use, see chapter 1.

4.2 Connecting the disposable stimulation probes

The Dr. Langer Medical GmbH disposable stimulation probes are fitted with DIN 42802-1 compliant connectors. The bipolar probes have two connectors (black and white), while the monopolar probes have only one connector (black). When using monopolar disposable stimulation probes, an extra subdermal needle electrode is required for stimulation as reference electrode (Art. No.: 42-0057).

4.2.1 Connection to the AVALANCHE® neuromonitors and the TWISTER® light neurostimulator

You need one of the following adapters to connect the disposable stimulation probes to the AVALANCHE® neuromonitors and to the TWISTER® light neurostimulators of Dr. Langer Medical GmbH:

Adapter	Article No.	Connector
for AVALANCHE® headbox, AVALANCHE® PLUS I/O box and TWISTER® light	41-0048, 41-0080	Connector green/grey
for AVALANCHE® PLUS I/O box	41-0087	Connector grey/grey

The adapter has two DIN 42802-1 sockets for connection of the disposable stimulation probes and for connection to the stimulation sockets of the AVALANCHE® neuromonitors and the TWISTER® light neurostimulator.


Figure 1: Adapter (left: Art. No.: 41-0048, right: Art. No. 41-0087)
for connection of the disposable stimulation probes to equipment of Dr. Langer Medical GmbH

Bipolar disposable stimulation probes

When using the bipolar disposable stimulation probes (Art. No. 40-1304, 40-1307, 40-1308), the two connectors of the probe (black and white) are connected to the sockets of the adapter cable. No reference electrode is required for bipolar stimulation!

- Connect the appropriate adapter to the green or grey stimulation socket of the AVALANCHE[®] neuromonitor or of the TWISTER[®] light neurostimulator.
- Connect the two plugs of the bipolar disposable stimulation probe to the sockets of the respective adapter! Be sure to observe the colour codes of the connectors (see Figure 1)!

Monopolar disposable stimulation probes

When using the monopolar disposable stimulation probes (Art. No. 40-1303, 40-1306), an extra disposable needle electrode is required for stimulation as reference electrode (Art. No. 42-9712); this needle electrode is to be inserted into the tissue near the stimulation site. The needle electrode is inserted at an angle and should subsequently be secured to prevent it from slipping out inadvertently. The needle electrode also has a DIN 42802-1 connector (white) and it is connected to the respective adapter cable.

- Connect the appropriate adapter to the green or grey stimulation socket of the AVALANCHE[®] neuromonitor or of the TWISTER[®] light neurostimulator.
- Connect the black plug of the monopolar disposable stimulation probe to the socket of the adapter cable marked with a black dot (see Figure 1)!
- Stick the sterile subdermal needle electrode (Art. No. 42-9712) into the surrounding tissue near the stimulation site and secure the needle electrode!
- Connect the white plug of the needle electrode to the white socket of the adapter (see Figure 1)!

NOTICE:

Do not autoclave the adapter cables (Art. No. 41-0048, 41-0080, 41-0087); autoclaving may damage the adapter! Keep the adapters outside the sterile area!

4.3 Connection to the TWISTER[®] MM neurostimulator

To connect the disposable stimulation probes to the TWISTER[®] MM neurostimulator of Dr. Langer Medical GmbH, you need one of the following stimulation cables:

Stimulation cable	Article No.	Connector
for TWISTER [®] MM output 1 - 4	41-0073	Connector blue/black
for TWISTER [®] MM output 5 - 6	41-0083	Connector black/black

The adapters have two DIN 42802-1 sockets for connection of the plugs of the disposable stimulation probes and one plug for connection to the stimulation sockets of the TWISTER[®] MM neurostimulator.


Figure 2: Stimulation cable (Art. No.: 41-0073)
for connection of the disposable stimulation probes to equipment of Dr. Langer Medical GmbH

Bipolar disposable stimulation probes

When using the bipolar disposable stimulation probes (Art. No. 40-1304, 40-1307, 40-1308), both connectors of the probe (black and white) are plugged into the sockets of the stimulation cable. No reference electrode is required for bipolar stimulation!

- Connect the appropriate stimulation cable to the blue or black stimulation socket of the TWISTER[®] MM neurostimulator.
- Connect the two plugs of the bipolar disposable stimulation probe to the sockets of the respective stimulation cable! Be sure to observe the colour codes of the connectors (see Figure 2)

Monopolar disposable stimulation probes

When using the monopolar disposable stimulation probes (Art. No. 40-1303, 40-1306), an extra disposable needle electrode is required for stimulation as reference electrode (Art. No. 42-9712); this needle electrode is to be inserted into the tissue near the stimulation site. The needle electrode is inserted at an angle and should subsequently be secured to prevent it from slipping out inadvertently. The needle electrode also has a DIN 42802-1 connector (white) and it is connected to the respective stimulation cable.

- Connect the appropriate stimulation cable to the blue or black stimulation socket of the TWISTER[®] MM neurostimulator.
- Connect the black plug of the monopolar disposable stimulation probe to the socket of the stimulation cable marked with a black dot (see Figure 2)
- Stick the sterile subdermal needle electrode (Art. No. 42-9712) into the surrounding tissue near the stimulation site and secure the needle electrode!
- Connect the white plug of the needle electrode to the white socket of the stimulation cable (see Figure 2)!

NOTICE:

Do not autoclave stimulation cables (Art. No.: 41-0073, 41-0083); autoclaving may damage the cables! Keep the stimulation cables outside the sterile area!

5 Handling

5.1 Unpacking

The Dr. Langer Medical GmbH disposable stimulation probes are delivered sterile. Remove the probe from its packaging under aseptic conditions.

5.2 Preparation and performance check

The disposable stimulation probes of Dr. Langer Medical GmbH consist of a handle with the probe tips and an integrated connection cable. Monopolar probes have one connector, whereas bipolar probes have two DIN 42802-1 compliant connectors.

Before using disposable stimulation probes, check the following points:

- The electrode tip of the probe must not be bent and must not exhibit any signs of mechanical deformation.
- The insulation of the connection cable, especially at the connectors, must not be damaged.
- The electrical contacts in the plugs must not be bent, the opposite contact openings in the sockets must be unobstructed.
- The connectors must mate easily. Insertion of the cable into the sockets must be possible without feeling any resistance.
- Before use, the plug on the connection cable of the disposable stimulation probe must be connected to the socket of the neuromonitor/neurostimulator in the non-sterile area.
- During surgery it must be ensured that neither the probe nor the connection cable can fall down into the non-sterile area.

5.3 Stimulation

Apply the electrode tip of the probe to the tissue that you intend to stimulate. With bipolar probes, both poles of the probe must be in contact with the tissue at the same time.

6 Troubleshooting

Description of the problem	Possible cause	Remedy
No or weak response of the muscle to stimulation	Poor contact between electrode and nerve / muscle	Remove blood or rinsing solution
	Low current intensity	Select higher current setting
	Muscle relaxant	Contact anaesthetist
	Defective probe	Replace stimulation probe

7 Disposal

 CAUTION	Disposal Defective and used devices may not be disposed of as unsorted municipal waste. • The device must be disposed of in compliance with the rules for the disposal of surgical waste valid in your facility and with local regulations and legislation for the disposal of surgical waste.
--	--

8 Technical specifications

Product-specific data	
Length of the connection cable	2.5 m
Connector	complies with DIN 42802-1
Standards	
Classification according to EU Directive 93/42/EEC	III

9 Warranty

The Terms and Conditions of Dr. Langer Medical GmbH or the conditions of the purchase agreement apply.

The warranty covers only the material.

Warranty claims are recognised only if the receiver/user proves that the following rules have been adhered to:

- Intended use of the disposable stimulation probes was observed.
- Proper handling.
- Instructions for transport, storage and application were observed.
- No unauthorised structural modifications were carried out.
- No unauthorised repairs were carried out.

Liability for consequential damage is excluded in any case.

10 Copyright and reproduction of the instructions for use

The copyright to the instructions for use remains with Dr. Langer Medical GmbH, Waldkirch. The instructions for use may not be copied without the approval of the manufacturer. All rights are reserved for these instructions for use, in particular the right to reproduction and distribution as well as translation. No part of these instructions for use may be reproduced in any form (e.g., photocopy, microfilm or other procedure), processed, copied or distributed using electronic systems without the prior written consent of Dr. Langer Medical GmbH.

The information contained in these instructions for use may be modified or expanded without prior notice and does not represent any obligation on the part of Dr. Langer Medical GmbH.

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Copyright © 2021
All rights reserved
Dr. Langer Medical GmbH Şartları ve Koşulları veya satın alma sözleşmesinin koşulları geçerlidir.