Dr. Langer Medical GmbH

Am Bruckwald 26
79183 Waldkirch
GERMANY

Fon: +49 (0) 76 81 / 47 45 4-0
Fax: +49 (0) 76 81 / 47 45 4-29

www.medical-langer.de
info@medical-langer.de

DR. LANGER
M E D I C A L

SAXOPHONELEKTRODE® Gebrauchsanweisung	DEUTSCH	2
SAXOPHONE ELECTRODE Instructions for Use	ENGLISH	8
ÉLECTRODE SAXOPHONE Notice d'utilisation	FRANÇAIS	14
ELETTRODO A SASSOFONO Istruzioni per l'uso	ITALIANO	21
ELECTRODO SAXOFÓN Instrucciones de uso	ESPAÑOL	27
ELÉTRODO TIPO SAXOFONE Instruções de utilização	PORTUGUÊS	34
SAXOFOONELEKTRODE Gebruiksaanwijzingen	NEDERLANDS	40
ELEKTRODA „SAKSOFON“ Instrukcja obsługi	POLSKI	47
SAXOFONO VÁ ELEKTRODA Návod k použití	ČESKY	54
ΗΛΕΚΤΡΟΔΙΟ ΤΥΠΟΥ ΣΑΞΟΦΩΝΟΥ Οδηγίες χρήσης	ΕΛΛΗΝΙΚΑ	60
ЭЛЕКТРОД ТИПА «САКСОФОН» Инструкции по эксплуатации	РУССКИЙ	67
SAKSAFON ELEKTRODU Kullanım Kılavuzu	TÜRKÇE	74

Dr. Langer Medical GmbH
Am Bruckwald 26
79183 Waldkirch
GERMANY

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Instructions for use Art. No. 05-0071

SAXOPHONE ELECTRODE

Saxophone electrode for single use (Art. No. 42-0068)

CONTENTS, ENGLISH

1	Intended use	8
2	Safety notices and precautions.....	9
3	Symbols used on the packaging	10
4	Combination and connection of the saxophone electrode.....	11
5	Handling.....	12
6	Troubleshooting	13
7	Repairs and structural modifications.....	13
8	Disposal	13
9	Technical specifications	13
10	Warranty	13
11	Copyright and reproduction of the instructions for use	13

1 Intended use

1.1 Intended purpose

The single-use saxophone electrode by Dr. Langer Medical GmbH – referred to as saxophone electrode below – is intended for use in intraoperative neurological monitoring during surgical interventions for continuous stimulation of the vagus nerve. The saxophone electrode consists of a silicone carrier with integrated electrode surfaces, a 175-cm cable and two 1.5-mm DIN 42802-1 compliant connectors. The electrode is used in combination with the AVALANCHE® neuromonitor. The saxophone electrode is connected to the VANESSA headbox and to the neuromonitor with an adapter. The saxophone electrode, Art. No. 42-0068, is a single-use device and has been sterilized with ethylene oxide before shipment. Do NOT reprocess the device. After exposure of the vagus nerve, the saxophone electrode is placed in the carotid sheath between the two large vessels in the neck (V. jugularis interna, A. carotis communis) around the vagus nerve (for detailed handling and application instructions, refer to section "5"). This establishes an electrical contact with the vagus nerve below the two vessels. Via the electrode, the nerve can be stimulated directly. Based on the muscle action potential, which is elicited by the nerve stimulation, the function of the recurrent laryngeal nerve, branching off the vagus nerve, can be identified and verified. If the electrode is used as intended, its design prevents tissue injuries.

1.2 Indication

The saxophone electrode by Dr. Langer Medical GmbH is only intended for direct stimulation of the vagus nerve and, in combination with an EMG measuring device, for monitoring of the vagus nerve and the recurrent laryngeal nerve, branching off the vagus nerve. With this electrode, the entire nerve path of the recurrent laryngeal nerve can be monitored for integrity. The AVALANCHE® neuromonitors display the EMG signals of the muscle action potentials of the target muscles (vocal muscle) elicited by electrical stimulation, which can then be assessed by the physician/surgeon.

1.3 Contraindication

Do NOT use the saxophone electrode if the course of the vagus nerve is atypical and, in particular, if its position is ventral rather than dorsal with respect to the cervical vessels. The electrode could cause lesions to the nerve by exerting excessive pressure. The stimulation site must be kept clear of excessive amounts of blood or rinsing solution, because the electrical energy of the stimulation pulse would be dispersed across the liquid and, as a result, would not be sufficient for stimulation of the nerve. Furthermore, the tissue to stimulate should not be too dry, because a high electrical resistance at the site where the electrode is in contact with the tissue will lead to an inadequate stimulation current. The warning "Check vagus probe!" emitted by the AVALANCHE® neuromonitor may indicate this situation. The saxophone electrode is not suitable for minimally invasive surgery or surgical interventions with a very small incision. The saxophone electrode is not intended for direct contact with the heart or the central circulatory system. The saxophone electrode may not be used in pregnant patients, patients with active implants, such as a pacemaker or ICDs, or epileptic patients.

1.4 User profile

The saxophone electrode may only be used in patients for the indication in section 1.2 by physicians who meet the following requirements:

- They have the requisite training, knowledge, and experience.
- They have read and understood the instructions for use.
- They are experienced in handling sterile medical devices.
- They have an understanding of neurophysiology especially regarding intraoperative nerve and muscle stimulation and recording of muscle action potentials.
- They have knowledge of the application areas of intraoperative neuromonitoring and are aware of the contraindications.

- They know how to recognise and prevent false negative stimulation reactions during direct nerve stimulation,
- They are aware of the effects of muscle relaxants on EMG signals.

The saxophone electrode may only be prepared and passed on to the surgeon by OR staff who meet the following requirements:

- They have the requisite training, knowledge, and experience.
- They have read and understood the instructions for use.
- They are experienced in handling sterile medical devices.

1.5 Patient profile

In conjunction with neuromonitors and neurostimulators, the saxophone electrode by Dr. Langer Medical GmbH is applied during surgical interventions according to the indication in section 1.2. There are no known clinical limitations regarding age or gender of the patients.

1.6 Exclusions from the intended purpose

The saxophone electrode is not intended for therapeutic purposes. The saxophone electrode is not intended for use in neurodiagnostics or outside of the OR. The administration of muscle relaxants will partially or completely block the stimulus conduction from the motor nerve to the innervated muscle. Intraoperative neuromonitoring while using stimulators (identification and function monitoring of motor nerves, evoked potentials) cannot be fulfilled in this case.

1.7 Use, transport, storage and disposal

The saxophone electrode is intended for intraoperative application in the OR. For handling instructions, please refer to chapter 5 in these instructions for use. For detailed information about transport and storage of the saxophone electrode, please refer to the device labelling and to chapter 3 in these instructions for use. For information about disposal of the medical device, please refer to chapter 8. The saxophone electrode may not be disposed of as unsorted municipal waste. It must be disposed of in compliance with the regulations for the disposal of surgical waste applicable in your facility.

2 Safety notices and precautions

The requisite personal protective equipment must be provided by the operator. All existing safety devices must be checked regularly.

2.1 The meaning of safety notices

 DANGER	• indicates an imminent hazard that can result in death if it is not avoided. • indicates an imminent hazard that can result in devastation of the environment if it is not avoided.
 WARNING	• indicates a potential hazard that can result in serious bodily injuries if it is not avoided. • indicates a potential hazard that can affect the proper functioning of nearby equipment if it is not avoided. • indicates a potential hazard that can result in significant damage of the environment if it is not avoided.
 CAUTION	• indicates a potential hazard that can result in mild or moderate bodily injuries if it is not avoided. • indicates a potential hazard that can result in brief malfunctions of the device if it is not avoided.
ATTENTION	• indicates a potential hazard that can result in property damage if it is not avoided.

2.2 The meaning of notices

"Notice:"

- indicates manufacturer information which directly or indirectly concerns the safety of individuals or the protection of property. The information does not directly concern a hazard or a hazardous situation.
- indicates manufacturer information that is important or useful for the operation or maintenance of the device.

2.3 Knowledge of the instructions for use

Understanding the instructions for use is absolutely essential for correct operation of the product. Therefore, everyone who

- prepares
- works with
- and cleans/reprocesses

the product, must read the instructions for use. Observe the safety notices of each chapter.

2.4 Informal safety measures

These instructions for use must be readily accessible at the location where the devices are used. The safety notices and precautions must be kept in a legible condition.

2.5 Observing the safety notices

Working with medical device is generally associated with residual risks to patients, to the medical staff and to the environment. Risks cannot be entirely eliminated by design measures alone. Safety does not depend solely on the product. Safety depends largely on the training of medical staff and correct operation of the product. Everyone working with the device is expected to read, understand and observe the safety notices in this chapter.

2.6 Safety notices and precautions

The safety notices are divided into the following hazard categories: Defective device / defective accessories, application and clinical residual risks, anaesthesia.

2.6.1 Defective device, defective accessories

 CAUTION	Defective/modified device or connected accessories When using a defective/modified device and/or accessories, safe operation and handling cannot be guaranteed. The user and patient may suffer injuries. • Check the device and all accessories for damage before each use (mechanical damage, kinked cable, damaged insulation, bent or broken probe tips...). • A damaged device and/or accessory may not be used. Discard defective devices and accessories.
---	---

2.6.2 Application and clinical residual risks

 WARNING	Use with neuromonitors/neurostimulators not designed for this purpose Using the devices with neuromonitors/neurostimulators not designed for their intended purpose may cause high current densities at the electrode surfaces. There is a risk of skin burns and tissue lesions. • The device may only be employed in compliance with its intended use, see chapter 1. • Use the device only with appropriate neuromonitors/neurostimulators and observe their instructions for use.
 WARNING	Use of device from damaged pack/past its expiry date There is a risk of infection when a device from a damaged pack and/or after its expiry date is used. • Check the packaging for damage. • Check the expiry date.
 WARNING	Use of a reprocessed device The saxophone electrode is a disposable device. It is not intended to be reprocessed. The procedures employed for reprocessing the device with the purpose of reusing it, including its repeated sterilisation, may degrade its mechanical and biological characteristics to an unacceptable level and considerably increase the risk of contamination. There is a risk of infection, allergies, toxic reactions and other intolerances. • Do not reprocess the device. • After use, dispose of the device according to the instructions in chapter 8 of these instructions for use
 WARNING	Muscle does not respond although nerve/muscle is intact (false negative response to stimulation). The stimulation current actually applied to the nerve may be reduced by current division according to Ohm's law, particularly due to fluid accumulated in the area of the electrode-tissue contact area or when the stimulation electrode is too far away from the excitable structure. • Prevent fluid accumulation.
 WARNING	Stimulation not possible Dry tissue and/or dried tissue residue on the electrode surfaces causes a very high contact impedance between electrode and tissue. It is not possible to deliver a stimulation pulse. • Clean the electrode surfaces. • Moisten the tissue.
 WARNING	Simultaneous use of an electrosurgical unit and the device, or device in contact with HF electrode. Improper use of the electrosurgical unit (wrong connection of HF electrode and/or HF electrode in contact with electrode surfaces or connection cables of the saxophone electrode) may cause the transfer of high frequency currents to the patient. Risk of skin burns at reference electrode, lead electrode or stimulation electrode or other conductive parts. • Do not allow the device to come into contact with the HF electrode. • Observe the instructions for use of the electrosurgical unit.
 WARNING	Continuous vagus stimulation affecting organs controlled by autonomic nervous system In rare situations, the continuous vagus stimulation may affect autonomic functions regulating the heart rate, for example, and lead to cardiac arrhythmias. • Contact the anaesthetist.
 WARNING	Use of muscle relaxants The administration of muscle relaxants will partially or completely block the stimulus conduction from the motor nerve to the innervated muscle. In this case, the device cannot fulfil its intended purpose, i.e., the localisation of nerves by means of electrical stimulation. • Do not use the device if muscle relaxants are administered. • Contact the anaesthetist.

3 Symbols used on the packaging

	Identifies the manufacturer of the medical device according to the EU Directive 93/42/EEC.		Identifies a double sterile barrier system.
	Article number. Indicates the manufacturer's order number which allows the medical device to be identified.		Identifies a double sterile barrier system with protective packaging outside.

	Lot code. Indicates the manufacturer's lot code, which allows the batch or lot to be identified.		Identifies a disposable or single-patient use/one-time-use medical device.
	Indicates the date of manufacture.		Identifies a medical device which does not contain any natural rubber latex.
	Indicates the expiry date after which the use of the medical device is no longer permitted.		Identifies a medical device.
	Indicates a medical device which was sterilised with ethylene oxide.		Unique Device Identification
	Attention. Indicates that the user is required to study the instructions for use, particularly with respect to safety information such as warnings and precautions that, for various reasons, cannot be placed directly on the medical device.	 (01)0425180720 0643 (11)210100 (17)230100 (10)P000000-01	GS1 Datamatrix 2D code 2D code as legible text: (01) Global Trade Item Number (GTIN) (11) Date of manufacture (17) Expiry date (10) LOT/batch code
	Observe instructions for use. Indicates that the user is required to refer to the instructions for use.		Symbol indicates compliance with the applicable EU Directives. 0297: ID number of the notified body (DQS Medizinprodukte GmbH, Frankfurt am Main, Germany)
	Indicates a medical device which should not be used if its packaging is damaged or open.		
	Identifies a medical device which must be protected from light sources.		Identifies a medical device which must be protected from humidity.

4 Combination and connection of the saxophone electrode

4.1 Combination with accessories and other medical devices

Accessories and combinable medical devices by Dr. Langer Medical GmbH are recommended for use with the saxophone electrode. All information about connections and technical specifications of accessories and medical devices provided in this document refer to products from Dr. Langer Medical GmbH.

 WARNING	Use with neuromonitors/neurostimulators not designed for this purpose Using the devices with neuromonitors/neurostimulators not designed for their intended purpose may cause high current densities at the electrode surfaces. There is a risk of skin burns and tissue lesions. • The device may only be employed in compliance with its intended use, see chapter 1. • Use the device only with appropriate neuromonitors/neurostimulators and observe their instructions for use.
---	---

4.2 Connection to the AVALANCHE® neuromonitors

The Dr. Langer Medical GmbH saxophone electrode is fitted with DIN 42802-1 compliant connectors (white and black). You need one of the following adapters to connect the saxophone electrode to the AVALANCHE® neuromonitors of Dr. Langer Medical GmbH:

Adapter	Article No.	Connector
for AVALANCHE® VANESSA headbox	41-0048	Connector green/grey
	41-0080	

The adapters have two DIN 42802-1 sockets for connection of the plugs of the saxophone electrode and one plug for connection to the stimulation sockets of the AVALANCHE® neuromonitors.

- Connect the green/grey connector of the adapter to the socket labelled "Output 2 Vagus", which is the right socket on the VANESSA headbox (pre-amplifier) of the AVALANCHE® neuromonitor (see Figure 1).
- Plug the DIN connector of the saxophone electrode into the sockets of the adapter. Be sure to observe the colour codes of the connectors (see Figure 1).


Figure 1: Connecting the saxophone electrode (Art. No. 42-0068) to the adapter (Art. No. 41-0080) and the VANESSA headbox

5 Handling

5.1 Unpacking

The Dr. Langer Medical GmbH saxophone electrode is delivered in a sterile packaging system (EO sterilized). Remove the electrode from its packaging under aseptic conditions. Do NOT reprocess the device.

5.2 Preparation and inspection

The saxophone electrode consists of a silicone carrier with integrated electrode surfaces, a 175-cm cable and two 1.5-mm DIN 42802-1 compliant connectors.

Before using the saxophone electrode, check the following points:

- The saxophone electrode and the integrated electrode surfaces must not show any signs of damage (e.g., loose connections).
- The insulation of the connection cable, especially at the connectors, must not be damaged.
- The electrical contacts in the plugs must not be bent, the corresponding contact openings in the sockets must be unobstructed.
- The connectors must mate easily. Insertion of the plugs into the sockets must be possible without feeling any resistance.
- During surgery the saxophone electrode should be accommodated to prevent the electrode and the connection cable from falling down into the non-sterile area.

5.3 Application of the saxophone electrode

The saxophone electrode is placed between the two large vessels in the neck (common carotid artery, internal jugular vein), in the carotid sheath. The electrode must be placed in a position cranial to the point where the recurrent laryngeal nerve branches off the vagus nerve so as to allow the recurrent nerve to be stimulated indirectly. The electrode must be applied to the ipsilateral vagus nerve.

- Expose the carotid sheath (Figure 2) above the point where the recurrent laryngeal nerve branches off the vagus nerve.
- Open the carotid sheath with a small incision of about 1 to 2 cm to allow the saxophone electrode to be placed between the cervical vessels.
- To begin with, use the handheld probe to stimulate the vagus nerve. The saxophone electrode can be inserted as soon as an evoked signal is displayed on the screen of the AVALANCHE® neuromonitors.

NOTICE:

If an initial vagal signal cannot be elicited, the reasons may be:

- prolonged effects of relaxants
- electrode not properly connected to the device
- incorrect position of tube
- non-recurrent laryngeal nerve on the right – rare (approx. 0.5 %)¹
- Using a suitable forceps applicator, introduce the saxophone electrode dorsally to the vagus nerve. Then pull the electrode around the vagus nerve to position the nerve directly on the concave surface of the electrode (Figure 2 and Figure 3). To achieve stable and reproducible neuromonitoring, the saxophone electrode must encircle the nerve as completely as possible.


Figure 2: Schematic illustration of the saxophone electrode correctly placed in the carotid sheath (schematic cross-section)

Its special design allows the electrode to be securely held in place and to prevent unwanted stimulation of the surrounding structures. The electrical pulses are directly delivered to the vagus nerve.


Figure 3: Schematic representation a) application of the saxophone electrode by hooking it on the nerve; b) saxophone electrode placed correctly on the nerve; c) the saxophone electrode can be detached by pulling in the direction of the connection cable

¹ Brauckhoff M, Walls G, Brauckhoff K, Thanh PN, Thomusch O, Dralle H: Identification of the non recurrent inferior laryngeal nerve using intraoperative neurostimulation. *Langenbeck's Arch Surg* 386 (2000) 482-487

5.4 Displaced saxophone electrode

The special design of the saxophone electrode ensures that it is reliably held in place.

CAUTION

Incorrect placement of the device

The electrode may shift due to manipulations at the cable or surgical site. The AVALANCHE® neuromonitors monitor the current that flows between the poles of the saxophone electrode. If no current can be detected, the monitors emit an audio signal. To prevent any disturbing noise, the audio signal of the monitor is automatically suppressed after the alert.

- Review and correct the placement of the electrode when you hear the audio signal.

5.5 Removing the saxophone electrode

At the end of the surgical intervention, the saxophone electrode can be removed from the carotid sheath by gently pulling on the connection cable (Figure 3).

6 Troubleshooting

Description of the problem	Possible cause	Remedy
No or weak response of the muscle to stimulation	Poor contact between electrode and nerve / muscle due to blood or rinsing solution	Remove blood or rinsing solution
	Displaced saxophone electrode	Reposition the electrode
	Low current intensity	Select higher current setting
	Muscle relaxant	Contact anaesthetist
	Defective saxophone electrode	Use another electrode

7 Repairs and structural modifications

Unauthorised repairs and structural modifications are not allowed.

8 Disposal

CAUTION

Disposal

Defective and used devices may not be disposed of as unsorted municipal waste.

- The device must be disposed of in compliance with the rules for the disposal of surgical waste valid in your facility and with local regulations and legislation for the disposal of surgical waste.

9 Technical specifications

Product-specific data	
Length of the connection cable	175 cm
Connector	complies with DIN 42802-1
Standards	
Classification according to EU Directive 93/42/EEC	IIa

10 Warranty

The Terms and Conditions of Dr. Langer Medical GmbH or the conditions of the purchase agreement apply. The warranty covers only the material. Warranty claims are recognised only if the receiver/user proves that the following rules have been adhered to:

- Intended use of the disposable stimulation probes was observed.
- Proper handling.
- Instructions for transport, storage and application were observed.
- No unauthorised structural modifications were carried out.
- No unauthorised repairs were carried out.

Liability for consequential damage is excluded in any case.

11 Copyright and reproduction of the instructions for use

The copyright to the instructions for use remains with Dr. Langer Medical GmbH, Waldkirch. The instructions for use may not be copied without the approval of the manufacturer. All rights are reserved for these instructions for use, in particular the right to reproduction and distribution as well as translation. No part of these instructions for use may be reproduced in any form (e.g., photocopy, microfilm or other procedure) or processed, copied or distributed using electronic systems without the prior written consent of Dr. Langer Medical GmbH. The information contained in these instructions for use may be modified or expanded without prior notice and does not represent any obligation on the part of Dr. Langer Medical GmbH.

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Copyright © 2021
All rights reserved