Dr. Langer Medical GmbH

Am Bruckwald 26
79183 Waldkirch
GERMANY

Fon: +49 (0) 76 81 / 47 45 4-0
Fax: +49 (0) 76 81 / 47 45 4-29

www.medical-langer.de
info@medical-langer.de

DR. LANGER

TUBUSKLEBEELEKTRODE Gebrauchsanweisung	DEUTSCH	2
ADHESIVE TUBE ELECTRODE User Manual	ENGLISH	13
ÉLECTRODE TUBE ADHÉSIVE Instructions d'utilisation	FRANÇAIS	25
ELETTRODO ADESIVO PER TUBO ENDOTRACHEALE Istruzioni per l'uso	ITALIANO	37
ÉLECTRODE TUBE ADHÉSIVE Manual de usuario	ESPAÑOL	49
YAPIŞKAN TÜP ELEKTRODU Kullanım Kılavuzu	TÜRK	61

Standard-Ausführung: (Art.-Nr.: 42-9800)

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CE 0051

EN ISO 13485

Manufacturer: Spes Medica S.p.A., Via Buccari 21, 16153 Genova, Italy

User Manual Art. No.: 05-0061

All rights are reserved for this user manual, in particular the right to reproduction and distribution as well as translation. No part of this user manual may be reproduced in any form (e.g., photocopy, microfilm or other procedure) or processed, copied or distributed using electronic systems without the prior written consent of Dr. Langer Medical GmbH.

The information contained in this user manual may be modified or expanded without prior notice and does not represent any obligation on the part of Dr. Langer Medical GmbH.

Similar images.

Printed by Dr. Langer Medical GmbH

Printed in Germany

Copyright © Dr. Langer Medical GmbH, Waldkirch 2023

Contents

1	Intended Use	15
2	Indication	15
3	Contraindication	16
4	Safety Notices and Precautions	16
4.1	Observing the Safety Notices	16
4.2	Safety Notices and Precautions	16
5	Handling	18
5.1	Packaging of the adhesive tube electrode	18
5.2	Different adhesive tube electrode models	18
5.3	Attaching the adhesive tube electrode	18
5.4	Intubation	19
5.5	Connecting the adhesive tube electrode to the neuromonitor AVALANCHE®	19
5.6	Tap test	21
6	Disposal	22
7	Symbols and Explanations	22
8	Technical Specifications	24
9	Manufacturer and Distributor	24

1 Intended Use

The Dr. Langer Medical GmbH adhesive tube electrode is a disposable product intended for use in intra-operative neuromonitoring. During surgical interventions, it is used to record the muscle action potentials at the vocal cord muscles.

The adhesive tube electrode consists of the adhesive electrode coated with silver/silver chloride and of the attached connection cables. In addition, adhesive electrodes with push-button connectors (Art. No. 42-0016) and the corresponding push-button cable (Art. No. 41-0024) from Dr. Langer Medical GmbH are required. The connectors of the connection cables are DIN 42802-1 compliant.

The adhesive tube electrode is supplied sterile. The adhesive electrodes with push-button connectors are non-sterile and must not be attached to areas within the aseptic surgical site. The adhesive tube electrode as well as the adhesive electrode with push-button connector are disposable accessories. It is not permitted to reprocess and reuse them.

The adhesive tube electrode is suitable for limited contact duration of up to 24 hours (according to DIN EN ISO 10993-1).

The adhesive tube electrode must not be used to deliver energy (e.g., for electrical stimulation)!

NOTICE:

Only skilled persons trained in the use of the adhesive tube electrode are allowed to use it. Persons employing the adhesive tube electrode should be familiar with neuromonitoring functions, with nerve stimulation and nerve monitoring procedures. The operating facility undertakes to allow only persons who have been trained in the use of the adhesive tube electrode and who have read and understood the Safety Notices and Precautions chapter to work with this electrode.

2 Indication

The Dr. Langer Medical GmbH adhesive tube electrode is used for the acquisition of muscle action potentials at the vocal cord muscles with the aim of identifying the vagal nerve, the recurrent laryngeal nerve, or the superior laryngeal nerve and verifying their integrity. The adhesive tube electrode is also suitable for application in intraoperative neuromonitoring during craniotomy, carotid endarterectomy, skull base surgery, and anterior cervical fusion to protect the nerves at risk of being injured.

3 Contraindication

No contraindications for the use of adhesive tube electrodes are known to date.

4 Safety Notices and Precautions

Understanding the user manual is absolutely essential for the correct use of the product. Therefore, all persons preparing and using the product are requested to read the user manual carefully.

Please pay particular attention to the safety notices in each chapter.

4.1 Observing the Safety Notices

Working with medical products is generally associated with residual risks to the patient, the medical staff and the environment. Not all risks can be completely eliminated by design measures implemented by the manufacturer.

Safety does not depend solely on the equipment and on accessories. To a large extent, safety depends on the training of medical personnel, on the correct operation of the equipment and of the accessories and on the application in line with the intended use.

All persons working directly or indirectly with the product are required to read, understand, and observe the safety notices.

4.2 Safety Notices and Precautions

The use of muscle-relaxant anaesthetics influences the occurrence, amplitude and latency of muscle action potentials and may even suppress them altogether. Therefore, prior to the operation, the surgeon is required to thoroughly discuss the strategy and in particular the objectives of the surgical intervention, with the anaesthesia personnel.

All cables connected to the neuromonitor and to the accessories must be checked for a tight fit. Do not use cables with damaged or poorly fitting connectors.

Do not use connection cables if the insulation is damaged.

All devices and accessories must be checked before use and must not be applied if found to be damaged.

	Adhesive tube electrodes must not be used to deliver energy (e.g., for electrical stimulation).
	The adhesive tube electrode as well as the adhesive electrode with push-button connector are disposable products, which may not be reprocessed. Reprocessing of electrodes may damage or modify the materials used. Proper functioning is no longer guaranteed.
	The adhesive tube electrode as well as the adhesive electrode with push-button connector are disposable products, which may not be used after the expiry date. Modifications of the materials used may occur. Proper functioning is no longer guaranteed.
	The adhesive tube electrode is supplied sterile and must not be used if the packaging is damaged.
	Only to be used by the clinical specialist.
	It is necessary for the user and/or patient to report any serious incident occurring, related to the device, to the manufacturer and the competent authority of the member state where the user and/or patient is established.
	The environment of use is a hospital, specifically an operating room.
	The device is intended for both adult and pediatric patient populations. The device does not have limitations of use related to age (from birth onward), gender, or weight. The clinical specialist will identify the correct size for the clinical need.
	There are no contraindications or warnings for pregnant or lactating women.

It is essential to observe all special precautions and warnings concerning the use of Dr. Langer Medical GmbH adhesive tube electrodes with AVALANCHE® XT/SI neuromonitors! Furthermore, the instructions for use/safety notices/precautions of the stimulation devices must be observed.

Dr. Langer Medical GmbH adhesive tube electrodes are intended for use with neuromonitoring equipment. For use in conjunction with the adhesive tube electrodes, the AVALANCHE® neuromonitors from Dr. Langer Medical GmbH are recommended. All information concerning connections and technical specifications in this user manual refers to Dr. Langer Medical GmbH devices. If the Dr. Langer Medical GmbH adhesive tube electrode is used in conjunction with other neuromonitors, the user must observe the instructions for use of the neuromonitor concerned and, in particular, the intended use, the warning, and precautions. In all cases the neuromonitor used must be a medical device tested and approved in accordance with EN 60601!

The unauthorized replacement of the Dr. Langer Medical GmbH connectors with other connectors is strictly prohibited.

5 Handling

The reliable acquisition of the muscle action potentials by means of the adhesive tube electrode depends on the correct application of the electrode on the tube as well as on the exact positioning of the tube in the patient's glottis.

Make sure that no pull is exerted on the cable and do not secure the cables, not even by winding them around the tube.

5.1 Packaging of the adhesive tube electrode

The Dr. Langer Medical GmbH adhesive tube electrode is delivered sterile. Remove the electrode from its packaging under aseptic conditions.

5.2 Different adhesive tube electrode models

The following models of the Dr. Langer Medical GmbH adhesive tube electrode are available:

Model	Suitable tube size (inside diameter)	Article No.
standard model	di = 7,5 - 9 mm	42-9800

5.3 Attaching the adhesive tube electrode

- Check the inside diameter of the endotracheal tube and select the suitable adhesive tube electrode according to the table in section 5.2.
- First of all, remove the larger of the two application papers (marked with the number 1) from the back of the adhesive tube electrode.
- Align the electrode in the middle of the dorsal (outer) curvature of the tube – centre line of the tube = centre line of the electrode (see Figure 1 a).
- A distance of approx. 1 cm between cuff and electrode must be observed (see Figure 1 a).
- Use the excess length of the application paper to determine the correct distance. The distance is correct when the edge of the paper touches the cuff.
- Press the large adhesive surface completely onto the tube (see Figure 1 b).
- Remove the small application paper (marked with the number 2) and press the small electrode surface onto the tube so that it covers the transparent part of the larger electrode surface (see Figure 1 c).
- The two silver-coloured electrode surfaces must not touch each other or overlap.
- A correctly placed adhesive tube electrode should have as few creases as possible and be located in the middle of the tube.


Figure 1 Application of the adhesive tube electrode

- a) Lining up the centre lines of the tube and of the electrode and determining the distance of approx. 1 cm to the cuff by means of the protruding application paper
- b) Applying the larger section of the electrode
- c) Applying the smaller section of the electrode

5.4 Intubation

Attention: If lubricant is used for introduction of the tube, please apply the electrode first.

- Introduce the tube under visual control until the vocal cords contact the electrode surfaces centrally on both sides.
- Verify the electrode position and check whether the electrode surfaces are visible above the vocal cords.
- Place the tube in the centre line of the pharynx behind the tongue and tape it to the upper lip. The tube's centre line should be located directly in the gap between the front upper incisors.
- Secure the respirator tubes so that they will not twist the tube.

5.5 Connecting the adhesive tube electrode to the tube electrode to the AVALANCHE® neuromonitor


Plug the blue cable of the adhesive tube electrode into socket „A“ of the EMG acquisition cable (Art. No. 41-0010) and the red cable into socket „R“ as shown in Figure 2 and Figure 3.

Figure 2:
EMG cable with connected plugs


- Connect the adhesive electrode with push-button connector (Art. No. 42-0016), which is attached to the patient's forehead, cheek or neck, to the Dr. Langer Medical GmbH push-button cable (Art. No. 41-0024) and plug the green connector of the cable into the $\perp$ socket (see Figure 2 and Figure 3) of the EMG acquisition cable (Art. No. 41-0010).

Figure 3: Connecting the push-button electrode to the push-button cable and the adhesive tube electrode to the EMG cable

- Connect the EMG acquisition cable (Art. No. 41-0010) to yellow socket of the headbox (preamplifier) labelled „Channel 1...2“ or „Channel 1...4“ of the neuromonitor AVALANCHE® (see Figure 4).

Figure 4: Connecting the EMG cable to the ...

NOTICE:

The adhesive tube electrodes can also be connected to the neuromonitor AVALANCHE® with the following EMG acquisition cables and the corresponding headboxes.

Table 1 Connection options for the adhesive tube electrode

EMG signal acquisition cable Article No.	SI Headbox (preamplifier) Article No.	XT Headbox (preamplifier) Article No.
41-0010, 41-0012, 41-0068	21-0022-21-0028	21-0006-21-0008, 21-0013
41-0014, 41-0069	21-0022, 21-0023, 21-0025, 21-0026, 21-0027, 21-0028	21-0007, 21-0008
41-0016, 41-0070	21-0023, 21-0027	21-0008

5.6 Tap test


In addition to verifying the position visually, the electrode location can be assessed by means of the „tap test“ (see Figure 5).

A strong tap on the patient's throat to the right or left of the larynx should evoke peaks (deflections) in the free running EMG on the screen of the neuromonitor (see Figure 5).


Figure 5 Tap test

a) Tapping on the patient's throat

b) EMG signal without taps (baseline)

c) EMG signal with peaks (deflections) during the tap test

6 Disposal

Used Dr. Langer Medical GmbH adhesive tube electrodes must be disposed of in compliance with the regulations for the disposal of surgical waste applicable in your facility! It is not allowed to reuse or reprocess adhesive tube electrodes.

7 Symbols and Explanations

The user manual must be readily accessible at the location where the devices are used. The hazard and safety notices must be kept in a legible condition.

The following terms and symbols are used in the user manual and on the product packaging:

 WARNING	Indicates a potential hazard which can result in death or severe bodily injuries if it is not avoided.
 CAUTION	Indicates a potential hazard which can result in mild or moderate bodily injuries if it is not avoided.
CAUTION	Indicates a potential hazard which can result in property damage if it is not avoided.
NOTICE	Indicates manufacturer information which directly or indirectly concerns the safety of individuals or the protection of property. The information does not directly concern a hazard or a hazardous situation.
	Ethylene Oxide Sterilization with single sterile barrier system
	Ethylene Oxide Sterilization with single sterile barrier system with protective packaging outside.
	Do not reuse
	Keep dry

	Protect from sunlight
	Do not use if packaging is damaged
	Symbol indicates Complies with the Regulation (EU) 2017/745: ID of the Notified Body
	Observe operating instructions!
	Indicates that the user is required to study the user manual, particularly with respect to safety information such as warnings and precautions which, for various reasons, cannot be placed directly on the medical device.
	Use by
	Manufacturer
	Medical Device
	Unique Device Identification
	Swiss Authorized Representative

8 Technical Specifications

Product-specific data

Length of the connection cable

75 cm

Ambient conditions for storage and transport

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Copyright © 2023 Dr. Langer Medical GmbH
Copyright © Dr. Langer Medical GmbH, Waldkirch (Alemania) 2023