Natus Disposable Snap Electrodes

Instructions for Use:

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Mode d'emploi :

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EU Representative / Importer
Natus Manufacturing Limited
IDA Business Park
Gort, Co. Galway, Ireland

Arazy Group Switzerland GmbH
Bruderholzallee 53
4059 Basel
Switzerland
swiss.ar@arazygroup.com

UK Responsible Person
Natus Nicolet UK Ltd
Baynards Green Trading Estate
Prospect House
Oxfordshire
Bicester
OX27 7SG
England, United Kingdom

Rx only

Associated product part numbers: 1410500

031529 Rev D 12/2022

English

Description:

Snap electrodes are used in conjunction with snap leads and amplifiers for monitoring electromyogenic and physiologic data signals commonly recorded from chest, extremities and faces of patients.

Intended Use:

The Disposable Snap Electrodes are intended to be used with compatible 1.5 mm touchproof button snap cables for transferring surface biopotential signals during Sleep Studies, EEG, EMG, and recording evoked potential signals from a cutaneous location. For use on adults.

Intended User and Target Patient Group:

The product is for use only under the direction of trained medical professionals in acquiring electrophysiologic signals.

Clinical Benefits:

Facilitates obtaining the recording of an EEG/PSG study to detect any irregularities indicative of various brain or sleep disorders.

Contraindications and Side Effects:

There are no known contraindications or side effects for procedures performed with Natus Disposable Snap Electrodes. See Warnings or Precautions below.

Operating Instructions:

- Open the pouch just before use.
- Prepare the electrode site prior to application making sure that the site is clean and dry.
- Connect the appropriate snap lead to the snap connector on the electrode.
- Remove the electrode from its backing.
- Apply to the patient making sure to place even pressure across the entire electrode area.
- Electrodes can be used to record or stimulate.
- Re-seal pouch after each use.
- Once recording is complete, detach and discard the disposable snap electrode.

Understanding Cautions Statements:

CAUTION

Refers to a hazardous situation that could result in minor or moderate injury or material damage if not avoided.

- Information on how the hazardous situation is avoided.

Warnings and Precautions:

CAUTION

Device dropped or damaged in transit/use could lead to loss of function or delayed diagnosis.

- Inspect the device prior to each use and do not use if damaged.

Electrode reused on another patient leads to cross infection or loss of performance.

- Do not reuse or sterilize.

Data cannot be recorded due to improper placement leading to delay in procedure.

- Assure proper connections and check signal quality before using.

Snap electrodes left in place for prolonged periods of time may be associated with skin redness.

- Do not use the electrodes for prolonged periods of time.

Warnings and Precautions Continued:

CAUTION

Repositioning of electrode can degrade the signal quality.

- Electrode should be replaced if the adhesive no longer provides a uniform connection to the skin surface.

Modifications to the device can affect function and performance.

- Do not modify this equipment without authorization of the manufacturer.

Environmental Specifications:

Operating Conditions:

- Temperature: +10°C (+50°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 70 kPa to 106 kPa

Storage Conditions:

- Temperature: 0°C (+32°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 50 kPa to 106 kPa

Disposal Instructions:

Following use, dispose of adhesive electrodes with infectious waste.

Disclaimer:

Natus Medical Incorporated DBA Excel-Tech Ltd. (Xltek) is not responsible for injury, infection or other damage resulting from the use of this product.


Any serious incident that has occurred in relation to the device should be reported to Natus Medical Incorporated DBA Excel-Tech Ltd. (Xltek) and the competent authority of the Member State in which the user and/or patient is established.

Refer to the Natus website for an electronic copy of this document.

Glossary of Symbols:

Symbol	Standard Reference	Standard Title of Symbol	Standard Title as per Referenced Standard	Explanation
	-	-	An indication of Medical device	This product is a medical device.
Rx only	21 CFR Part 801.109(b)(1)	Labeling-Prescription devices.	Prescription only	Indicates the product is authorized for sale by or on the order of a licensed healthcare practitioner.
	ISO 15223-1 Symbol 5.4.5 (Reference Annex B for the general prohibition symbol)	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Not made with Natural Rubber Latex	Indicates a medical device that is not made with natural rubber latex.
	ISO 15223-1 Symbol 5.1.1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Manufacturer	Indicates the medical device manufacturer.
	ISO 15223-1 Symbol 5.1.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Authorized representative in the European Community	Indicates the Authorized representative in the European Community.
	ISO 15223-1 Symbol 5.1.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Symbol 5.1.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Batch or Lot code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Symbol 5.1.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Symbol 5.4.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	ISO 60601-1 Table D.1 #11	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Operating instructions	
	ISO 60601-1 Table D.2 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Follow instructions for use	Refer to instruction manual/ Booklet. NOTE on ME EQUIPMENT “Follow instructions for use”
	ISO 15223-1 Symbol 5.4.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Caution: Read all warnings and precautions in instructions for use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 60601-1 Table D.1 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.		

Glossary of Symbols Continued:

Symbol	Standard Reference	Standard Title of Symbol	Standard Title as per Referenced Standard	Explanation
	ISO 60601-1 Table D.2 #2	Medical electrical equipment — Part 1: General requirements for basic safety and Essential performance.	General warning sign	Indicates a hazard of potential personal injury to patient or operator.
	ISO 15223-1 Symbol 5.1.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Use-by date	Indicates the date after which the medical device is not to be used. Note: This symbol shall be accompanied by a date to indicate that the medical device should not be used after the end of the year, month or day shown.
	ISO 15223-1 Symbol 5.2.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1 Symbol 5.3.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Temperature limit	Indicates the (storage) temperature limits to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.3.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Humidity limitation	Indicates the range of (storage) humidity to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.4.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Swiss Medical Device Ordinance (MedDO)	Swiss Medical Device (MedDO).	Indicates the Authorized Representative in Switzerland	Indicates the Authorized Representative in Switzerland.
	UKCA Medical Device Regulation (SI 2022 No 618, as amended) (UK MDR 2002)	UKCA Medical Device Regulation.	UKCA Mark	In compliance with the United Kingdom technical conformity.
	EU Medical Device Regulations 2017/745	EU Medical Device Regulation.	CE marking	Signifies European technical conformity.

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