Tremor Sensor

Instructions for Use

Natus Manufacturing Limited
IDA Business Park
Gort, Co. Galway, Ireland
Tel: +353 (0)91 647400
Website: natus.com

Arazy Group Switzerland GmbH
Bruderholzallee 53
4059 Basel
Switzerland
swiss.ar@arazygroup.com

UK Responsible Person
Optima Medical Ltd (Neuro).,
30a Alan Turing Road
The Surrey Research Park
Guildford Surrey
GU2 7AA United Kingdom

Rx only

Associated Product Part Number:

9013L0241

Copyright © 2022-2026 by Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.

035978_05 04/2026

English

Description:

The Tremor Sensor is placed on the area of a tremor, typically the patient's hands, fingers, or arms, to provide movement data to the electrodiagnostic system. The data is used to analyze the characteristics of tremors associated with movement disorders.

Intended Use:

The Tremor Sensor is a device used to record the characteristics of tremor movements for the monitoring and evaluation of neurologic function.

Intended Users and Patient Target Group:

The Tremor Sensor is intended for use by skilled physicians and technologists trained in the specialty of Electrodiagnostic testing.

The target patient population is the pediatric and adult patient population requiring Electrodiagnostic testing.

Clinical Benefits:

The Tremor Sensor allows monitoring of a patient's tremor throughout an Electrodiagnostic test. This information can be useful criteria in interpretation of the test results. The Tremor Sensor is used at the discretion of the clinical provider.

Contraindications and Side Effects:

There are no known contraindications or side effects for procedures performed with the Tremor Sensor.

Operating Instructions:

Note: Inspect the device prior to each use and do not use if damaged.

1. Secure the Tremor Sensor to an appropriate location on the patient's skin using medical grade tape approved by your medical facility.
2. Plug the Tremor Sensor Connectors into the inputs of the EMG system amplifier.
3. Ensure that a Ground electrode is applied to the patient and connected to the amplifier for the EMG procedure.

Cleaning Instructions:

1. Wipe the sensor with a damp cloth and mild detergent.
2. Do not allow fluid to get into the connectors.
3. Gently wipe the sensor with a lint-free cloth and air dry.

Understanding Caution Statements:

	CAUTION
Refers to a hazardous situation that could result in minor or moderate injury or material damage if not avoided.	
• Information on how the hazardous situation is avoided.	

Cautions:

	CAUTION
Use of damaged cables may lead to loss of function and performance.	
• Inspect the cables before each use and check signal quality.	
Unauthorized modification could lead to loss of device function and performance.	
• Do not modify the device.	

Environmental Specifications:

Operating conditions:

- Temperature: +10°C to +35°C (+50°F to +95°F)
- Relative Humidity: 20% to 80%
- Atmospheric Pressure: 700 hPa to 1060 hPa

Storage/Transport conditions:

- Temperature: -25°C to +55°C (-13°F to +131°F)
- Relative Humidity: 10% to 95%
- Atmospheric Pressure: 700 hPa to 1060 hPa

Compliance Standards:

ISO 10993-1: 2018 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process

ETSI EN 300 019-2-1 Environmental Engineering (EE); Environmental conditions and environmental tests for telecommunications equipment; Part 2-1: Specification of environmental tests; Storage

ETSI EN 300 019-2-2 Environmental Engineering (EE); Environmental conditions and environmental tests for telecommunications equipment; Part 2-2: Specification of environmental tests; Transportation

ASTM D4169-16 Standard Practice for Performance Testing of Shipping Containers and Systems for Vibration

IEC 60601-1 Ed. 3.1 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance

Disposal Instructions:

Natus is committed to meeting the requirements of the European Union WEEE (Waste Electrical and Electronic Equipment) Regulations 2014. These regulations state that electrical and electronic waste must be separately collected for the proper treatment and recovery to ensure that WEEE is reused or recycled safely. In line with that commitment Natus may pass along the obligation for take back and recycling to the end user, unless other arrangements have been made. Please contact us for details on the collection

and recovery systems available to you in your region at natus.com

Electrical and electronic equipment (EEE) contains materials, components and substances that may be hazardous and present a risk to human health and the environment when WEEE is not handled correctly. Therefore, end users also have a role to play in ensuring that WEEE is reused and recycled safely. Users of electrical and electronic equipment must not discard WEEE together with other wastes. Users must use the municipal collection schemes or the producer/importers take-back obligation or licensed waste carriers to reduce adverse environmental impacts in connection with disposal of waste electrical and electronic equipment and to increase opportunities for reuse, recycling and recovery of waste electrical and electronic equipment.

Equipment marked with the below crossed-out wheeled bin is electrical and electronic equipment. The crossed-out wheeled bin symbol indicates that waste electrical and electronic equipment should not be discarded together with unseparated waste but must be collected separately.

Notice:

Any serious incident that has occurred in relation to the device should be reported to Natus Medical Incorporated and the competent authority of the Member State in which the user and/or patient is established.

Instructions to Access the eIFU:

An electronic copy of the Instructions for Use (eIFU) for this product is available to view and print from natus.com/support by navigating to the Product Knowledge Center, searching for the IFU part number, and choosing the version for your local language (as available). The files can be printed, saved, or searched using Adobe Reader which can be downloaded for free from Adobe Systems at adobe.com.

For delivery within seven business days, request a paper copy of the Instructions for Use in any of these languages by contacting Natus Technical Support or by contacting your local Natus Distributor/Representative. Contact information can be found in the eIFU. IFUs are subject to change without notice.

Glossary of Symbols:

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title	Explanation
	ISO 15223-1 Symbol 5.7.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	An indication of Medical Device	This product is a medical device.
Rx only	21 CFR Part 801.109(b)(1)	Labeling-Prescription devices	Prescription only	Indicates that the product is authorized for sale by or on the order of a licensed healthcare practitioner.
	2012/19/EU	Waste Electrical and Electronic Equipment (WEEE)	Disposal at end of operating life instructions	Indicates that electrical and electronic equipment waste should not be discarded together with unseparated waste but must be collected separately.
	ISO 15223-1 Symbol 5.1.1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Manufacturer	Indicates the medical device manufacturer.
	SO 15223-1 Symbol 5.1.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Symbol 5.1.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Symbol 5.1.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Symbol 5.4.3 Annex A #A.15	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Consult instructions for use	Indicates an instruction to consult an electronic instructions for use (eIFU).
	ISO 15223-1 Symbol 5.4.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	IEC 60601-1 Table D.1 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance		
	MDR 2017/745	EU Medical Device Regulation	CE marking	Signifies European technical conformity.
	ISO 15223-1 Symbol 5.3.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Temperature limit	Indicates the (storage) temperature limits to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.3.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Humidity limitation	Indicates the range of (storage) humidity to which the medical device can be safely exposed.

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title	Explanation
	ISO 15223-1 Symbol 5.3.9	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Atmospheric pressure limitation	Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.1.2	Swiss Medical Device Ordinance (MedDO)	Indicates the Authorized Representative in Switzerland	Indicates the Authorized Representative in Switzerland.
	UK MDR 2002	UKCA Medical Device Regulation (SI 2022 No 618, as amended)	UKCA Mark	In compliance with the United Kingdom technical conformity.
	ISO 15223-1 Symbol 5.1.11	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied	Country of Manufacture	Identifies the country of manufacture of products. The Date of Manufacture may appear next to this symbol.

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