

Precision Gold

DISPOSABLE EEG ELECTRODES

Disposable Gold 9mm Cup EEG Electrode, 50" (1.5m), 10pack

REF 13-10100-6010 PN SE-251-P QTY 1

LOT 0106-H10K EXP 2025-07-16



(01)00817377022852(11)250716(10)0106-H10K

DyMedix
DIAGNOSTICS



DyMedix Diagnostics, Inc
8985 Rice Creek Parkway
Shoreview, MN 55126 USA
888-212-1100
www.dymedix.com

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 - Ste 123
2585 AM The Hague
The Netherlands

EC REP

Instructions:

Inspect electrodes prior to use. Prepare the scalp using conventional technique. Fill the cup with a conductive paste recommended: Ten20™. After use, dispose of electrodes. These electrodes are not to be reused. Do not apply electrodes to broken or damaged skin.

Not made with Natural Latex Rubber.



Consult Instructions for Use electronic labeling at www.dymedix.com/downloads



Single Use Only

Rx Only (U.S.)

86-51800-0100 (B)