



# STERILE BLUNTED NEEDLE



## Description

Disposable sterile blunted needle that can be used to prepare the patient's skin and to insert the conductive gel into the prewire electrode housing during the EEG test.

It is recommended to discard the needles after each use according to hospital infection control procedures.

STERILE BLUNTED NEEDLE  
SPECIFIC DEVELOPED FOR  
EEG APPLICATIONS

UMDNS 12737

## Materials

Cannula	Stainless Steel
Needle cover	PP
Connection	PP
Packaging	Medical Paper/PP

## Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Sterilization method	Ethylene oxide (EtO)
Latex	No
Phthalates	No

## Specifications

Generic code	DISPXXXXXXXX
Dimensions	From 15 to 25 mm
Needle diameter	1,60 mm (16G)
Connector	Standard Luer Lock

## Environment

PVC-Free	✓
RoHS II	✓
REACH	✓