



Instructions For Use & Product Information

Disposable Adhesive Surface Electrodes



Caution/warnings:

- In the USA Federal Law restricts this device to sale by or on the order of a physician or licensed practitioner.
- Only to be used by trained and skilled personnel, under responsibility of a physician or licensed practitioner.
- If irritation develops, discontinue use and consult your clinician.
- The electrodes are non sterile.
- The size, shape and type of electrode may affect the safety and effectiveness of electrical stimulation and recording.
- The electrical performance characteristics of electrodes may affect the safety and effectiveness of electrical stimulation and recording.
- Using stimulation electrodes that are too small or incorrectly applied could result in discomfort or skin burns.
- The user should contact the manufacturer of the electrical stimulator or recording device if the user does not know if the electrode can be used with the stimulation or recording device.
- Do not share electrodes with other persons because of the risk of adverse skin reactions and disease transmissions.
- Storage temperature between 15 °C(59 °F) and 35 °C(95 °F).

Instructions for application:


- Apply only to unbroken skin. Thoroughly cleanse application site to remove lotion, dirt, etc. Dry site completely. Clip hair if necessary. **DO NOT SHAVE OR ABRASE AREA!**
- Remove electrodes from package. Remove electrodes from the liner by lifting the edge of the electrode.
- Apply electrodes to treatment site by firmly pressing down centre of pad first and then smoothing down to electrode edges.

Instructions for removal:

- Turn off device.
- Disconnect lead wire from device.
- Remove electrode from skin.
- Remove possible residue of the gel on the skin with some water and a normal soap.




USA
Rx Only



Instructions For Use & Product Information

Disposable Adhesive Surface Electrodes



Technomed
medical accessories

Caution/warnings:



- In the USA Federal Law restricts this device to sale by or on the order of a physician or licensed practitioner.
- Only to be used by trained and skilled personnel, under responsibility of a physician or licensed practitioner.
- If irritation develops, discontinue use and consult your clinician.
- The electrodes are non sterile.
- The size, shape and type of electrode may affect the safety and effectiveness of electrical stimulation and recording.
- The electrical performance characteristics of electrodes may affect the safety and effectiveness of electrical stimulation and recording.
- Using stimulation electrodes that are too small or incorrectly applied could result in discomfort or skin burns.
- The user should contact the manufacturer of the electrical stimulator or recording device if the user does not know if the electrode can be used with the stimulation or recording device.
- Do not share electrodes with other persons because of the risk of adverse skin reactions and disease transmissions.
- Storage temperature between 15 °C(59 °F) and 35 °C(95 °F).

Instructions for application:

- Apply only to unbroken skin. Thoroughly cleanse application site to remove lotion, dirt, etc. Dry site completely. Clip hair if necessary. **DO NOT SHAVE OR ABRASE AREA!**
- Remove electrodes from package. Remove electrodes from the liner by lifting the edge of the electrode.
- Apply electrodes to treatment site by firmly pressing down centre of pad first and then smoothing down to electrode edges.

Instructions for removal:

- Turn off device.
- Disconnect lead wire from device.
- ~~Remove electrode from skin.~~
- Remove possible residue of the gel on the skin with some water and a normal soap.



Instructions

Print : One sided black
Paper : Offset coat pressed 70 grs.
Finished : Folded
Packed : Manageable